


## Instructions for Use and Technical Description



### ELEGANZA 4

### Positionable Bed for Intensive Care

*version with scales and without scales*



D9U001GE4-0101

Version: 14

Publication Date: 2024-09

## **Manufacturer:**

L I N E T spol. s r.o.  
Želevčice 5  
274 01 Slaný

Tel.: +420 312 576 111  
Fax: +420 312 522 668

E-mail: [info@linet.cz](mailto:info@linet.cz)  
<http://www.linet.com>  
Service department: [service@linetgroup.com](mailto:service@linetgroup.com)

## **Authorized Representative in Great Britain:**

LINET UK Ltd  
11 Brunel Way  
Segensworth East  
Fareham  
Hampshire  
PO15 5TX  
United Kingdom

## **Authorized Representative in Switzerland:**

Bigla Care AG  
Bernstrasse 3  
CH-3421 Lyssach  
Switzerland

## **Authorized Representative in Malaysia:**

Emergo Malaysia Sdn. Bhd.  
Level 16, 1 Sentral, Jalan Stesen Sentral 5  
KL Sentral  
50470 Kuala Lumpur  
Malaysia

## **Authorized Representative in Peru:**

Drogería Emergo Peru S.R.L.  
Calle Las Orquídeas Nro. 585, Int. 1301  
San Isidro, Lima  
Director Técnico Renato Delgado Rivera  
RUC: 205 52 75 65 35  
Peru

## **Authorized Representative in Philippines:**

Biomedica Healthcare Inc.  
Unit 2103 City land 10, Tower 1  
Ayala Corner, Dela Costa Street  
Makati City, Metro Manila  
Philippines

**Eleganza 4**

Positionable bed for intensive care

Author: L I N E T, s.r.o.

Related links: [www.linnet.com](http://www.linnet.com)

**D9U001GE4-0101**

Version: 14

Publication Date: 2024-09

Copyright © L I N E T, s.r.o., 2024

Translation © L I N E T, 2024

All rights reserved.

All trademarks and brands are the property of the appropriate owners. The manufacturer reserves the right to changes in the contents of this instructions for use that relate to the product's technical regulations. It is for this reason that the contents of this instructions for use may indicate differences from the current manufacture of the product. Reproduction, also excerpts, only with prior permission of the publisher. Subject to changes due to technical developments. All technical data are rated data and are subject to construction and manufacturing tolerances.

## Table of Contents

|                                                                            |           |                                                          |           |
|----------------------------------------------------------------------------|-----------|----------------------------------------------------------|-----------|
| <b>1 Symbols and Definitions</b>                                           | <b>6</b>  | 11.5.1 Backrest                                          | 53        |
| 1.1 Warning Notices                                                        | 6         | 11.5.2 Thighrest                                         | 54        |
| 1.1.1 Types of Warning Notices                                             | 6         | 11.5.3 Mechanical Calfrest Adjustment                    | 56        |
| 1.1.2 Structure of Warning Notices                                         | 6         | 11.5.4 Bed Height                                        | 57        |
| 1.2 Instructions                                                           | 6         | 11.5.5 Autocontour                                       | 59        |
| 1.3 Lists                                                                  | 6         | 11.5.6 Examination Position                              | 60        |
| 1.4 Symbols on the Package                                                 | 7         | 11.5.7 Emergency Trendelenburg Position                  | 62        |
| 1.5 Symbols and Labels on the Bed                                          | 8         | 11.5.8 Anti-Trendelenburg Tilt and Trendelenburg Tilt    | 63        |
| 1.6 Acoustic signalisation                                                 | 13        | 11.5.9 CPR Position                                      | 65        |
| 1.7 Visual signalisation                                                   | 14        | 11.5.10 Mechanical Bed Extension                         | 66        |
| 1.7.1 NIGHT LIGHT                                                          | 14        | 11.5.11 Cardiac Chair Position                           | 69        |
| 1.7.2 Mains Power LED (Nurse Control Panel, Attendant Control Panel)       | 14        | 11.5.12 Mobilisation Position                            | 70        |
| 1.7.3 Accumulator indicator (Nurse Control Panel, Attendant Control Panel) | 14        | 11.5.13 Ergoframe                                        | 70        |
| 1.7.4 Lock LED (Nurse Control Panel, Attendant Control Panel)              | 15        | <b>12 Scales Control (only version with scales)</b>      | <b>71</b> |
| 1.8 Definitions                                                            | 15        | 12.1 Preparation                                         | 71        |
| 1.9 Abbreviations                                                          | 16        | 12.2 Taring                                              | 71        |
| <b>2 Safety Instructions</b>                                               | <b>17</b> | 12.3 Displaying                                          | 71        |
| <b>3 Intended use</b>                                                      | <b>20</b> | 12.4 Hold Mode                                           | 72        |
| 3.1 User population                                                        | 20        | 12.5 Bed Overload                                        | 72        |
| 3.2 Contraindications                                                      | 20        | 12.6 Bed Underload                                       | 72        |
| 3.3 Operator                                                               | 20        | 12.7 Weighing in tilt                                    | 72        |
| <b>4 Product Description</b>                                               | <b>21</b> | 12.8 Zeroing Scales                                      | 72        |
| <b>5 Technical Specification</b>                                           | <b>22</b> | <b>13 Bed Exit Monitoring (only version with scales)</b> | <b>73</b> |
| 5.1 Identification of Applied Parts (Type B)                               | 22        | 13.1 Preparation                                         | 73        |
| 5.2 Scales                                                                 | 22        | 13.2 Activation of Bed Exit Monitoring                   | 73        |
| 5.3 Mechanical Specifications (Eleganza 4)                                 | 22        | 13.3 Monitored Zone                                      | 73        |
| 5.4 Environment Conditions                                                 | 23        | 13.4 BED EXIT ALARM                                      | 74        |
| 5.5 Electrical Specifications (Eleganza 4)                                 | 23        | 13.4.1 Alarm Volume                                      | 74        |
| 5.6 Electromagnetic Compatibility                                          | 24        | 13.5 PAUSE                                               | 74        |
| 5.6.1 Manufacturer instructions - electromagnetic emissions                | 24        | 13.6 Deactivation of Bed Exit Monitoring                 | 74        |
| 5.6.2 Manufacturer instructions - electromagnetic susceptibility           | 25        | <b>14 Equipment</b>                                      | <b>75</b> |
| <b>6 Use and Storage Conditions</b>                                        | <b>26</b> | 14.1 Accessory Rail with plastic hooks (optional)        | 75        |
| <b>7 Scope of Delivery and Bed Variants</b>                                | <b>27</b> | 14.2 Linen Shelf (standard)                              | 75        |
| 7.1 Delivery                                                               | 27        | 14.3 Mobi-Lift® (optional)                               | 76        |
| 7.2 Scope of Delivery                                                      | 27        | 14.4 Brake Signal (optional)                             | 77        |
| 7.3 Eleganza 4 Variants                                                    | 27        | 14.5 i-Brake® (optional)                                 | 77        |
| <b>8 Putting into Service</b>                                              | <b>28</b> | 14.6 Retractable Fifth Castor (optional)                 | 77        |
| 8.1 Accumulator Activation                                                 | 29        | 14.7 LINIS SafetyPort (optional)                         | 77        |
| 8.2 Head Board and Foot Board                                              | 30        | 14.8 i-Drive Power (optional)                            | 78        |
| 8.3 Mattress support platform                                              | 30        | 14.8.1 i-Drive Power System - Basic Description          | 78        |
| 8.4 Potential Equalisation                                                 | 31        | 14.8.2 Safety instruction for i-Drive Power              | 78        |
| 8.5 Before Use                                                             | 32        | 14.8.3 Specifications of Use                             | 78        |
| 8.6 Transport                                                              | 32        | 14.8.4 Manipulation                                      | 79        |
| 8.7 Firmware                                                               | 32        | 14.8.5 i-Drive Power Activation/Deactivation             | 80        |
| <b>9 Power Cable</b>                                                       | <b>33</b> | 14.8.6 Powered Drive                                     | 80        |
| <b>10 Accumulator</b>                                                      | <b>33</b> | 14.8.7 Braking                                           | 81        |
| 10.1 Replacing the accumulator                                             | 34        | 14.8.8 Free Drive                                        | 81        |
| 10.2 Removing the Bed from Use                                             | 35        | 14.8.9 Accumulator                                       | 81        |
| 10.3 Deactivating the Accumulator                                          | 35        | 14.8.10 Fault Signalization                              | 82        |
| <b>11 Manipulation</b>                                                     | <b>36</b> | 14.8.11 Light Indicators                                 | 82        |
| 11.1 Siderails                                                             | 36        | 14.8.12 Technical Specifications                         | 82        |
| 11.2 Castor Control                                                        | 41        | 14.8.13 Electrical specification                         | 83        |
| 11.3 CPR Backrest Release                                                  | 42        | 14.8.14 i-Drive Power Maintenance                        | 83        |
| 11.4 Control Elements                                                      | 43        | 14.9 Safestop (optional)                                 | 83        |
| 11.4.1 Nurse Control Panel (optional)                                      | 44        | 14.10 X-Ray Examination (optional)                       | 84        |
| 11.4.2 iBoard Basic (standard)                                             | 47        | 14.11 Undercarriage Cover                                | 85        |
| 11.4.3 Attendant Control Panel                                             | 49        | 14.12 USB Connector (optional)                           | 86        |
| 11.4.4 Handset (optional)                                                  | 51        | <b>15 Mattress</b>                                       | <b>87</b> |
| 11.4.5 Foot Switch - Bed Height (optional)                                 | 52        | 15.1 Installation of Passive Mattress                    | 87        |
| 11.4.6 Patient Control Panels                                              | 52        | 15.1.1 Straps with side release buckles                  | 87        |
| 11.5 Bed Positioning                                                       | 53        | 15.2 Installation of Active Mattress                     | 88        |
|                                                                            |           | 15.2.1 Straps with side release buckles                  | 88        |
|                                                                            |           | <b>16 Accessories</b>                                    | <b>90</b> |
|                                                                            |           | 16.1 Lifting Pole                                        | 92        |
|                                                                            |           | 16.2 Infusion Stand                                      | 93        |
|                                                                            |           | 16.3 Oxygen Bottle Holder                                | 94        |
|                                                                            |           | 16.4 Ventilation Circuit Holder                          | 95        |
|                                                                            |           | 16.5 Writing Shelf                                       | 96        |
|                                                                            |           | 16.6 Monitor Shelf                                       | 96        |

---

|                                                                            |            |
|----------------------------------------------------------------------------|------------|
| 16.7 Urinary Bag Holder .....                                              | 97         |
| 16.8 Protectors .....                                                      | 97         |
| 16.9 Safety Monitor .....                                                  | 100        |
| 16.10 Traction Frame E .....                                               | 101        |
| 16.11 Pushing Handles .....                                                | 103        |
| 16.12 Pushing Handles with Holder of the i-Drive Power Control Panel ..... | 104        |
| 16.13 IDock (size S and size M) .....                                      | 106        |
| <b>17 Cleaning/Disinfection .....</b>                                      | <b>108</b> |
| 17.1 Cleaning (Eleganza 4) .....                                           | 109        |
| 17.1.1 Daily Cleaning .....                                                | 109        |
| 17.1.2 Cleaning before Changing Patients .....                             | 109        |
| 17.1.3 Complete Cleaning and Disinfection .....                            | 110        |
| <b>18 Troubleshooting .....</b>                                            | <b>110</b> |
| <b>19 Maintenance .....</b>                                                | <b>111</b> |
| 19.1 Regular maintenance .....                                             | 111        |
| 19.2 Spare Parts .....                                                     | 111        |
| 19.3 Safety Technical Checks .....                                         | 111        |
| <b>20 Disposal .....</b>                                                   | <b>112</b> |
| 20.1 Environment Protection .....                                          | 112        |
| 20.2 Disposal .....                                                        | 112        |
| 20.2.1 Within Europe .....                                                 | 112        |
| 20.2.2 Outside Europe .....                                                | 112        |
| <b>21 Warranty .....</b>                                                   | <b>113</b> |
| <b>22 Standards and Regulations .....</b>                                  | <b>113</b> |

## 1 Symbols and Definitions

### 1.1 Warning Notices

#### 1.1.1 Types of Warning Notices

Warning notices are differentiated by the type of danger using the following key words:

- ▶ **CAUTION** warns about the risk of material damage.
- ▶ **WARNING** warns about the risk of physical injury.
- ▶ **DANGER** warns about the risk of fatal injury.

#### 1.1.2 Structure of Warning Notices



**SIGNAL WORDS!**

Type and source of danger!

- ▶ Measures to avoid the danger.

### 1.2 Instructions

Structure of instructions:




- ▶ Perform this step.
- Results, if necessary.

### 1.3 Lists




Structure of bulleted lists:

- List level 1
  - List level 2
    - List level 3


## 1.4 Symbols on the Package




|                                                                                     |                                                                       |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
|    | FRAGILE, HANDLE WITH CARE                                             |
|    | THIS WAY UP                                                           |
|    | KEEP DRY (PROTECT FROM HUMIDITY)                                      |
|    | PAPER RECYCLING SYMBOL                                                |
|   | OVERSEAS PACKAGE: STACKING LIMIT BY NUMBER (4 PACKAGES FOR STORAGE)   |
|  | OVERSEAS PACKAGE: STACKING LIMIT BY NUMBER (2 PACKAGES FOR TRANSPORT) |
|  | DO NOT USE HAND TRUCK HERE                                            |
|  | DO NOT STACK DURING STORAGE                                           |



## 1.5 Symbols and Labels on the Bed

|                                                                                                 |                                                             |
|-------------------------------------------------------------------------------------------------|-------------------------------------------------------------|
|                | READ INSTRUCTIONS FOR USE                                   |
|                | GO BUTTON (PRESS TO ACTIVATE CONTROL ELEMENT)               |
|                | STOP BUTTON (PRESS TO INTERRUPT BED POSITIONING)            |
|  <b>250 kg</b> | SAFE WORKING LOAD                                           |
|                | DESIGNATION OF HOSPITAL BED FOR ADULTS                      |
|              | USE MATTRESS RECOMMENDED BY MANUFACTURER                    |
|              | DO NOT PUT ANY OBJECTS ON UNDERCARRIAGE                     |
|              | BED EXTENSION                                               |
|              | WARNING AGAINST CRUSHING OR TRAPPING                        |
|              | JACK FOR ATTACHMENT OF CONDUCTOR FOR POTENTIAL EQUALISATION |
|              | CPR LEVER                                                   |
|              | WARNING                                                     |

|                                                                                     |                                                                                   |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
|    | WARNING                                                                           |
|    | MAXIMUM CALFREST LOAD                                                             |
|    | APPLIED PARTS TYPE <b>B</b>                                                       |
|    | THERMAL PROTECTION FOR TRANSFORMER                                                |
|    | SAFETY ISOLATING TRANSFORMER, GENERAL                                             |
|   | ONLY SUITABLE FOR INDOOR USE                                                      |
|  | MAXIMUM WEIGHT OF PATIENT                                                         |
|  | WEIGHT OF BED                                                                     |
|  | CE MARKING FOR ELEGANZA 4 WITHOUT SCALES                                          |
|  | CE MARKING FOR ELEGANZA 4 WITH SCALES                                             |
|  | WEEE SYMBOL (RECYCLE AS ELECTRONIC WASTE,<br>DO NOT PUT INTO THE HOUSEHOLD WASTE) |
|  | MANUFACTURER                                                                      |

|                                                                                               |                                                                                                                                                                                   |
|-----------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|              | DO NOT POLLUTE THE ENVIRONMENT                                                                                                                                                    |
|              | RECYCLING SYMBOL                                                                                                                                                                  |
|              | DO NOT OPEN                                                                                                                                                                       |
|              | ON                                                                                                                                                                                |
|             | OFF                                                                                                                                                                               |
|            | BATCH NUMBER (ACCESSORIES)                                                                                                                                                        |
| <br>Pb     | THIS COMPONENT CONTAINS LEAD - DO NOT PUT INTO THE HOUSEHOLD WASTE<br>(Eleganza 4 bed is equipped with lead accumulator. )                                                        |
| <br>Li-ion | THIS COMPONENT CONTAINS LITHIUM - DO NOT PUT INTO THE HOUSEHOLD WASTE<br>(If Eleganza 4 is equipped with Integration Module (SafetyMonitor System), it contains lithium battery.) |
|            | MEDICAL DEVICE (compatible with Medical Device Regulation)                                                                                                                        |

|                                                                                     |                                                            |
|-------------------------------------------------------------------------------------|------------------------------------------------------------|
|    | MANUFACTURING DATE                                         |
|    | REFERENCE NUMBER (PRODUCT TYPE DEPENDING ON CONFIGURATION) |
|    | SERIAL NUMBER                                              |
|    | DO NOT WASH (PROTECTOR COVER)                              |
|   | DO NOT BLEACH (PROTECTOR COVER)                            |
|  | DO NOT TUMBLE DRY (PROTECTOR COVER)                        |
|  | DO NOT WRING (PROTECTOR COVER)                             |
|  | DO NOT IRON (PROTECTOR COVER)                              |
|  | PROFESSIONAL CHEMICAL CLEANING (MODERATE PROCEDURE)        |
|  | DO NOT USE PHENOL (PROTECTOR COVER)                        |

|                                                                                     |                                                                                                                               |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
|    | DRY NATURALLY (PROTECTOR COVER)                                                                                               |
|    | RINSE WITH WATER (PROTECTOR COVER)                                                                                            |
|    | DISINFECT USING SOLUTION CONTAINING LESS THAN 1000 ppm OF CHLORINE (REFER TO INSTRUCTIONS FOR USE) (PROTECTOR COVER)          |
|    | HANDWASH WITH DETERGENT (INITIAL TEMPERATURE OF HOT WATER SHOULD NOT EXCEED 50°C) (PROTECTOR COVER)                           |
|  | COVER MATERIALS ARE FIRE RESISTANT TO BS7175, SOURCE 0, 1 AND 5 (PROTECTOR COVER)                                             |
|  | MASS OF MOBILE HOSPITAL BED (MASS OF EMPTY BED + SAFE WORKING LOAD)                                                           |
|  | UK CONFORMITY ASSESSED (UKCA) MARKING (ONLY ELEGANZA 4 WITHOUT SCALES NORMATIVELY HARMONIZED FOR GREAT BRITAIN ECONOMIC AREA) |
|  | AUTHORIZED REPRESENTATIVE IN GREAT BRITAIN                                                                                    |
|  | AUTHORIZED REPRESENTATIVE IN SWITZERLAND                                                                                      |

## 1.6 Acoustic signalisation

| SOUND                                                                  | MEANING                                                                                    |
|------------------------------------------------------------------------|--------------------------------------------------------------------------------------------|
| CONTINUOUS SOUND                                                       | overheating                                                                                |
|                                                                        | accumulator overcurrent                                                                    |
|                                                                        | scales overload (only version with scales)                                                 |
|                                                                        | actuator overload                                                                          |
| REPEATED BEEP: 0,6s sound / 2,6s silence                               | STOP error (all STOP buttons are disabled)                                                 |
| MELODY: 3 beeps, pause, 2 beeps, longer pause, 3 beeps, pause, 2 beeps | Bed Exit Alarm (only version with scales)                                                  |
| BEEP lasting 0,3s                                                      | confirmation                                                                               |
|                                                                        | stopping or locked function                                                                |
|                                                                        | optionally: transition from tilt (Trendelenburg, Antitrendelenburg) to horizontal position |
| BEEP lasting 0,5s                                                      | lowering to the lowest position                                                            |
|                                                                        | start of service mode or end of service mode                                               |
|                                                                        | keyboard error (positioning blocked)                                                       |
| BEEP lasting 3s                                                        | system error                                                                               |
| BEEP lasting 5s                                                        | scale module disconnected (only version with scales)                                       |
| REPEATED BEEP during 3 minutes: 1,1s sound / 1,1s silence              | Brake Signal (only version with Brake Signal)                                              |

## 1.7 Visual signalisation

### 1.7.1 NIGHT LIGHT

Bed illumination helps the nursing staff as well as the patient to orientate. The lowered intensity of lightning is set up after turning the bed on. The night light is turned off during accumulator operation.

The bed is equipped with three-phase illumination:

1. Lowered intensity of illumination
2. Full intensity of illumination
3. Illumination is turned off

#### After pressing any button:

► The bed illumination will light up at full intensity for 10 minutes.

After 10 minutes the bed illumination will be lowered.

After disconnection of the bed from the mains illumination lights up for few seconds.

#### Turning off bed illumination:

► Disconnect bed from the mains.

After disconnection of the bed from the mains illumination lights up for few seconds.

### 1.7.2 Mains Power LED (Nurse Control Panel, Attendant Control Panel)

| MAINS POWER LED<br> | MEANING                                                               |
|------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
| lit                                                                                                  | connected to the mains                                                |
| flashing: 0,6s lit / 0,6s not lit                                                                    | keyboard error (flashing inverted to Lock LED)<br>error (first fault) |
| flashing: 0,1s lit / 0,1s not lit                                                                    | service mode                                                          |
| not lit                                                                                              | disconnected from the mains power<br>transformer switching error      |

### 1.7.3 Accumulator indicator (Nurse Control Panel, Attendant Control Panel)

| ACCUMULATOR INDICATOR<br>  | MEANING                            |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|
| lit                                                                                                                                                                                              | accumulator disconnected or faulty |
| flashing: 1,6s lit / 0,2s not lit                                                                                                                                                                | accumulator deeply discharged      |
| flashing: 0,1s lit / 0,1s not lit                                                                                                                                                                | accumulator discharged             |
| flashing: 0,2s lit / 1,6s not lit                                                                                                                                                                | accumulator is charging            |
| not lit                                                                                                                                                                                          | accumulator charged                |

## 1.7.4 Lock LED (Nurse Control Panel, Attendant Control Panel)

| <br>LOCK LED | VISUAL<br>SIGNALISATION |  | lit    | 2 flashing modes in the case of<br>the control element with single<br>lock button:<br>1) flashing: 0,2s lit / 0,9s not lit<br>2) flashing: 0,1s lit / 0,1s not lit | flashing:<br>0,6s lit / 0,6s not lit |                   |                   | not lit  |
|-----------------------------------------------------------------------------------------------|-------------------------|--|--------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|-------------------|-------------------|----------|
|                                                                                               |                         |  |        |                                                                                                                                                                    |                                      |                   |                   |          |
| Thighrest Lock LED                                                                            |                         |  | locked | 1) lock settings - not selected lock<br>2) lock settings - selected lock                                                                                           | lock<br>error                        | keyboard<br>error | motion<br>blocked | unlocked |
| Backrest Lock LED                                                                             |                         |  | locked | 1) lock settings - not selected lock<br>2) lock settings - selected lock                                                                                           | lock<br>error                        | keyboard<br>error | motion<br>blocked | unlocked |
| Bed Height, Trendelenburg and An-<br>titrendelenburg Tilt Lock LED                            |                         |  | locked | 1) lock settings - not selected lock<br>2) lock settings - selected lock                                                                                           | lock<br>error                        | keyboard<br>error | motion<br>blocked | unlocked |
| Foot Switch Lock<br>LED                                                                       |                         |  | locked | 1) lock settings - not selected lock<br>2) lock settings - selected lock                                                                                           | lock<br>error                        | keyboard<br>error | motion<br>blocked | unlocked |

## 1.8 Definitions

|                             |                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Basic Bed Configuration     | the pricelist model configuration, not including a mattress                                                                                                                                                                                                                                                                                                                      |
| Bed Weight                  | The value depends on the product configuration, accessories or customer adjustments.                                                                                                                                                                                                                                                                                             |
| Clearance of Undercarriage  | the height from the floor to the lowest point of the undercarriage between the castors, for the manipulation of accessories under a braked bed in the standard position                                                                                                                                                                                                          |
| Duty Cycle                  | cycle of operation of the motor: time of activity/time of rest                                                                                                                                                                                                                                                                                                                   |
| Ergoframe                   | Ergoframe is the kinematic system of mattress support platform Adjustment whose effect is the elimination of pressure on the patient's abdomen and pelvic area and frictional forces on the patient's back and legs.                                                                                                                                                             |
| Maximum Patient Weight      | Maximum Patient Weight depends on the application environment according to IEC 60601-2-52. For application environment 1 (intensive/critical care) and 2 (acute care) reduce Safe working Load by 65 kg. For application environment 3 (long-term care) and 5 (ambulatory care) reduce Safe working Load by 35 kg.                                                               |
| Safe Working Load           | the highest allowable load on the bed (patient, mattress, accessories and the load supported by those accessories)                                                                                                                                                                                                                                                               |
| Siderail Height             | the height of the upper crossbar or the edges of the siderails (not the highest point of the siderail controls) from the patient surface                                                                                                                                                                                                                                         |
| Standard Bed Position       | <ul style="list-style-type: none"> <li>- The height of the patient surface with regard to the floor is 400 mm</li> <li>- The mattress support platform, including the individual parts, has to be in a horizontal (level - 0°) position.</li> <li>- The siderails are always locked in the upper position.</li> <li>- The basic position of the integrated extension.</li> </ul> |
| Adult                       | Patient having a physical size equal to or more than 146 cm, a mass equal to or more than 40 kg and a body mass index (BMI) equal to or more than 17 (according to IEC 60601-2-52).                                                                                                                                                                                              |
| Mass of mobile hospital bed | Sum of empty bed mass and Safe Working Load.                                                                                                                                                                                                                                                                                                                                     |

## 1.9 Abbreviations

|                   |                                                            |
|-------------------|------------------------------------------------------------|
| <b>AC ( ~ )</b>   | Alternating Current                                        |
| <b>ACP</b>        | Attendant Control Panel                                    |
| <b>CE</b>         | European Conformity                                        |
| <b>CPR</b>        | Cardiopulmonary Resuscitation                              |
| <b>dBA</b>        | Sound Intensity Unit                                       |
| <b>DC ( --- )</b> | Direct Current                                             |
| <b>CUC</b>        | Configuration number                                       |
| <b>EMC</b>        | Electromagnetic Compatibility                              |
| <b>FET</b>        | Field-effect transistor                                    |
| <b>HF</b>         | High Frequency                                             |
| <b>HPL</b>        | High Pressure Laminate                                     |
| <b>HW</b>         | Hardware                                                   |
| <b>ICU</b>        | Intensive Care Unit                                        |
| <b>INT.</b>       | Duty Cycle                                                 |
| <b>IP</b>         | Ingress Protection                                         |
| <b>IV</b>         | Intravenous                                                |
| <b>LED</b>        | Light Emitting Diodes                                      |
| <b>ME</b>         | Medical Electrical (Equipment)                             |
| <b>ON</b>         | Activation                                                 |
| <b>OFF</b>        | Deactivation                                               |
| <b>ppm</b>        | parts per million, millionth (1000 ppm = 0,1%)             |
| <b>REF</b>        | Reference Number (product type depending on configuration) |
| <b>SCU</b>        | System Control Unit (active mattress)                      |
| <b>SN</b>         | Serial Number                                              |
| <b>SW</b>         | Software                                                   |
| <b>SWL</b>        | Safe Working Load                                          |
| <b>UDI</b>        | Unique Device Identification (for medical devices)         |
| <b>USB</b>        | Universal Serial Bus                                       |
| <b>WEEE</b>       | Waste Electrical and Electronic Equipment                  |

## 2 Safety Instructions



**WARNING!**

Eleganza 4 bed should be left in its lowest position when the patient is unattended in order to reduce risk of injury due to falls!



**WARNING!**

Siderails of Eleganza 4 should be located in the „up“ position to reduce the risk of the patient accidentally slipping or rolling off the mattress!



**WARNING!**

Incompatible siderails and mattresses can cause an entrapment hazard!



**WARNING!**

Inappropriate handling of the power supply cord, e. g. by kinking, shearing or other mechanical damages is hazardous!



**WARNING!**

When routing cables from other equipment in the Eleganza 4 bed avoid squeezing those between parts of the Eleganza 4 bed!



**WARNING!**

Eleganza 4 bed should not be used with bed hoists and bed lifts!



**WARNING!**

The bed is intended for adults.  
► Follow the chapter **Intended use**.



**WARNING!**

Incompatible mattresses can create hazards.



**WARNING!**

To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.



**WARNING!**

No modification of this equipment is allowed.



**WARNING!**

Do not modify this equipment without authorization of the manufacturer.



**WARNING!**

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.



**WARNING!**

An additional multiple socket-outlet or extension cord shall not be connected to the medical electrical system.



**WARNING!**

During specific investigations or treatments the significant risks of reciprocal interference posed by ME equipment may occur.



**WARNING!**

Staff expert assessment is needed to consider all individual cases of contraindications!



**WARNING!**

Certain bed positions are not suitable for specific diagnosis/medical conditions. Fowler position is not suitable for spinal cord injuries! Trendelenburg position is not suitable for patients with higher intracranial pressure!



**WARNING!**

Length adjustment of the bed must be proportional to the height of patient!  
Risk of trapping or squeezing because of patient's body constitution disproportionate to the size of mattress support platform!



**WARNING!**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established!



**WARNING!**

Only authorised and trained person using the tool is allowed to change fuses and power supplies!



**WARNING!**

This medical device is not intended for oxygen enriched environment!



**WARNING!**

This medical device is not intended for use with flammable substances!



**WARNING!**

This medical device is not portable medical electrical equipment!



**WARNING!**

Make sure the duty cycle (2 min ON/18 min OFF) is not exceeded during bed positioning!



**WARNING!**

Patient is allowed to use selected control elements only if hospital personnel had assessed that the patient's physical and psychological state is in accordance with use of them and only if the hospital personnel had trained the patient in accordance with the instructions for use!



**WARNING!**

Hospital personnel is allowed to use the weighing system (scales) for weighing patients only if they had been trained according to the instructions for use!

---

**Additional Instructions for Correct Use:**

- ▶ Follow the instructions for use carefully.
- ▶ Use the bed exclusively if it is in perfect working order.
- ▶ If necessary, check the bed functions daily or at each shift change.
- ▶ Ensure any user has read and understood the instructions for use completely before operating the product.
- ▶ Use the bed exclusively with the correct mains supply.
- ▶ Ensure that the bed is operated exclusively by qualified personnel who have been trained according to the instructions for use.
- ▶ Ensure that the patient (health permitting) has been informed about the operation of the bed and all applicable safety instructions.
- ▶ Move the bed exclusively on even, hard-surfaced floors.
- ▶ Contact service department of the manufacturer immediately to let it replace any damaged parts by the original spare parts.
- ▶ Ensure that maintenance and installation are performed exclusively by qualified personnel who have been trained by the manufacturer.
- ▶ During peak loads or unavoidable excess loads (CPR), adjust Mattress support platform to the lowest position.
- ▶ Ensure that only one adult patient lies on the bed at any time.
- ▶ To avoid injury or crushing, take extra caution when operating any moving parts of the bed.
- ▶ When using lifting poles or infusion stands, ensure that nothing will be damaged when you move or adjust the bed.
- ▶ Brake the castors when the bed is occupied.
- ▶ Keep the mattress support platform in the lowest position at any time when the healthcare personnel are not treating the patient in order to prevent the patient from falling or injuries.
- ▶ Enable or disable functions on Patient Control Panels using the Attendant Control Panel as appropriate for the patient's physical and mental state. Verify that the function is actually disabled.
- ▶ Ensure that siderails are operated exclusively by healthcare personnel.
- ▶ Never use the bed in areas where there is a hazard of explosion.
- ▶ Never handle the mains plug with wet hands.
- ▶ Disconnect the product from the mains exclusively by pulling the mains plug.
- ▶ When pulling the mains plug, always hold the plug, not the cable.
- ▶ Position the mains cable so that there are no loops or kinks in the cable; protect the cable from mechanical wear and tear.
- ▶ Improper handling of mains cable can cause an electric shock hazard and other serious injuries or damages.
- ▶ To prevent failures, use exclusively the manufacturer's original accessories and mattresses in perfect working order.
- ▶ Ensure that the stipulated safe working load is not exceeded.
- ▶ If the patient's condition could lead to an entrapment, leave the mattress support platform in the flat position whilst unattended.
- ▶ Adjust bed height when transporting the bed in order to facilitate overcoming possible obstacles.
- ▶ Do not exceed maximum load of 80 kg for bed extension.
- ▶ Do not exceed the maximum patient weight limit (see Mechanical Specifications).
- ▶ Do not hang anything on any cable.
- ▶ Select a suitable location for the placement of bed accessories and other objects to prevent involuntary activation of buttons or controls which may result in the adjustment of bed positioning.
- ▶ Do not use the bed when its parts have been removed (e.g. parts of mattress support platform) unless these parts are designed to be removed.
- ▶ Never place any accessories or handset on the siderails where keyboards are located.
- ▶ After each emergency situation always check if any of the controllers (in siderails, handset or ACP) is not involuntarily pressed by the bed accessories or by the mattress.
- ▶ The weighing system must be calibrated at regular intervals and in accordance with the metrological regulations of the relevant country. All testing and certification must be carried out by qualified personnel. The healthcare provider is responsible for ensuring the required testing frequency and testing procedure of the weighing system is carried out.
- ▶ To avoid unintended activation of moving parts during any use of the bed always check that none of the control elements of the bed is pressed by persons, mattress or other objects.
- ▶ Before setting the Extra Low position, ensure there is no risk of any parts of the bed colliding with servers, accessories or body parts.
- ▶ Close Linen Shelf before using the Anti-Trendelenburg Tilt or Cardiac Chair Position.

## 3 Intended use

**The intended use is the hospitalization of the patient in the intensive and acute care units, which includes above all the following aspects:**

- ▶ Adjustment of the specific positions needed for the preventive reasons, routine nursing, treatments, mobilization, physiotherapy, examinations, sleeping, and relaxation. These positions are further specified and described in the clinical evaluation of this device, together with their potential clinical outcomes and benefits.
- ▶ Providing the safe environment for the patient during all relevant procedures. The particular requirements on patient safety are the subject of the clinical evaluation, including evaluation of the risk/benefit ratio. The relevant safety issues are the part of the risk management file.
- ▶ Patient in-bed indoor transport out of the patient room.
- ▶ Providing the suitable working conditions for the caregivers to perform the routine and specific tasks during the patient hospitalization.
- ▶ Indicative measurement of the patient weight, used as supportive feature without direct diagnostic effect. It helps staff to assess the general patient status and apply the nutrition and medicaments (valid for the version of the beds with in-bed scales).

### 3.1 User population

- ▶ Adult patients (weight  $\geq 40$  kg, height  $\geq 146$  cm, BMI  $\geq 17$ ) in the intensive and acute care units (Application Environment 1 and 2 as in IEC 60601-2-52)
- ▶ Caregivers (nurses, doctors, technical personnel, transport personnel, cleaning personnel)


### 3.2 Contraindications

- ▶ The medical device is not intended for the pediatric patients use.
- ▶ Certain positions are not suitable for specific diagnoses/medical conditions (e.g. spinal cord injuries vs. Fowler position, higher ICP patients vs. Trendelenburg). Staff expert assessment / nursing consideration is needed in all individual case of contraindication.

### 3.3 Operator

- ▶ Caregiver
- ▶ Patient (based on individual patient status assessment by caregiver the patient can utilize dedicated device functions)

## 4 Product Description



*Fig. Bed Overview (Eleganza 4)*

1. Removable Head Board
2. Head Siderail
3. Mobi-Lift Handle - Bed Height
4. Patient Control Panel
5. Foot Siderail
6. Removable Foot Board
7. Corner Bumper
8. Foot Board Lock
9. Linen Holder
10. Castor Control Lever
11. Foot Switch - Bed Height
12. Four-part Mattress support platform with Ergoframe® System
13. Siderail Release Mechanism
14. Accessory Rail
15. CPR Lever – Backrest Release
16. Castor
17. Accessory Adaptor
18. Nurse Control Panel
19. iBoard Basic

## 5 Technical Specification

All technical data are rated data and are subject to construction and manufacturing tolerances.

### 5.1 Identification of Applied Parts (Type B)

All part of the bed (and accessories) the patient can reach are type B Applied Parts.

- Mattress support platform Frame, Covers and all Movable Parts
- Head and Foot End
- Siderails
- Mobi-Lift® Handles
- Handset
- Patient Control Panels

### 5.2 Scales

Accuracy of displayed weight values:

- 0,5 kg
- Scales Class III

### 5.3 Mechanical Specifications (Eleganza 4)

| Parameter                                                                        | Value              |
|----------------------------------------------------------------------------------|--------------------|
| External Dimensions in Standard Bed Position (length x width)                    | 217,5 cm x 100 cm  |
| Maximum Siderail Height above Mattress support platform                          | 45 cm              |
| Dimensions of Head Siderail (length x height)                                    | 51,1 cm x 46,8 cm  |
| Dimensions of Foot Siderail (length x height)                                    | 100,2 cm x 43,1 cm |
| Bed Extension                                                                    | 11 cm, 22 cm       |
| Maximum Dimensions of Mattress (length x width)                                  | 208 cm x 90 cm     |
| Maximum Mattress Height                                                          | 23 cm              |
| Clearance of Undercarriage                                                       | 15 cm              |
| Castor diameter                                                                  | 15 cm              |
| Minimum-Maximum Mattress support platform Height (above floor, without Mattress) | 41 cm - 79 cm      |
| Ergoframe (Backrest/Thighrest)                                                   | 7,4 cm / 4 cm      |
| Maximum Backrest Angle                                                           | 65°                |
| Maximum Thighrest Angle                                                          | 25°                |
| Maximum Calfrest Angle                                                           | 20°                |
| Angle between Calfrest and Thighrest                                             | 225°               |
| Trendelenburg Angle                                                              | 14°                |
| Anti-Trendelenburg Angle                                                         | 14°                |
| Bed Weight (depending on configuration)                                          | 195 Kg             |
| SWL (Bed Safe Working Load)                                                      | 250 Kg             |
| SWL (Lifting Pole Safe Working Load)                                             | 75 Kg              |
| Maximum Patient Weight                                                           | 185 Kg             |
| Mass Of Mobile Hospital Bed (mass of empty bed + safe working load)              | 476 Kg             |
| Sound Pressure Level                                                             | 57 dBA             |
| Minimum Sound Pressure Level of Bed Exit Alarm                                   | 60 dBA             |

## 5.4 Environment Conditions

| Use Conditions                   |                    |
|----------------------------------|--------------------|
| Ambient Temperature              | 10°C - 40°C        |
| Relative Humidity                | 30% - 75 %         |
| Atmospheric Pressure             | 795 hPa - 1060 hPa |
| Storage and Transport Conditions |                    |
| Ambient Temperature              | -20°C - 50°C       |
| Relative Humidity                | 20% - 90 %         |
| Atmospheric Pressure             | 795 hPa - 1060 hPa |

## 5.5 Electrical Specifications (Eleganza 4)

| Parameter                            | Value                                                                        |
|--------------------------------------|------------------------------------------------------------------------------|
| Input Voltage                        |                                                                              |
| Version 1                            | 230 V~, 50/60 Hz                                                             |
| Version 2                            | 100 V~, 50/60 Hz                                                             |
| Version 3                            | 110 V~, 50/60 Hz                                                             |
| Version 4                            | 120 V~, 50/60 Hz                                                             |
| Version 5                            | 127 V~, 50/60 Hz                                                             |
| Version 6                            | 110 - 127 V~, 60 Hz or 220 V~, 50/60 Hz                                      |
| Maximum Power Input                  | 370 VAC                                                                      |
| Ingress Protection (EN 60529)        | IP X4                                                                        |
| Protection Class                     | Class I                                                                      |
| Electrical Motor Duty Cycle          | 2 minutes ON /18 minutes OFF                                                 |
| Accumulator                          | Pb AKU 2 x 12 V / 1,2 Ah / Fuse 15 A                                         |
| Fuse (version without i-Drive Power) |                                                                              |
| Version 1                            | 2x T1.6A L 250 V for 230 V, 110-127 V or 220 V version                       |
| Version 2                            | 2x T3.15A L 250 V for 100 V, 110 V, 120 V, 127 V, 110-127 V or 220 V version |
| Fuse (version with i-Drive Power)    |                                                                              |
| Version 1                            | 2 x T2A L 250V for 230 V                                                     |
| Version 2                            | 2 x T4A L 250V for 100 V, 110 V, 120 V, 127 V                                |

**NOTE** Upon request, LINET ® can deliver hospital beds with electrical specifications that comply with regional standards (custom voltage, different mains plugs).

## 5.6 Electromagnetic Compatibility

Bed is intended for hospitals except for near active HF surgical equipment and the RF shielded room of a medical system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Bed has defined no essential performance.



### WARNING!

It is recommended to avoid the use of this device next to or in block with other device, because it could lead to improper operation. If such use is needed, this device and the other equipment should be under surveillance to verify proper operation.

### List of used cables:

- ▶ **Mains cable**, maximum length 6 m
- ▶ **Attendant Control Panel**, maximum length 3m
- ▶ **Handset**, maximum length 3m



### WARNING!

Use of the accessories, converters and other cables, than specified and provided by manufacturer of this bed could lead to increase of electromagnetic emission or lower the electromagnetic immunity of this bed and lead to improper operation.



### WARNING!

Mobile RF communication device (including end use devices like antenna cables and external antenna) should not be used closer than 30 cm (12 inches) from any part of this bed Eleganza 4, including cables specified by manufacturer. Otherwise this could lead to deterioration of functionality of this bed.



### WARNING!

Do not overload the bed (SWL), respect the duty cycle (INT.) and consider chapter 19 Maintenance in order to maintain the basic safety with regard to electromagnetic disturbances for the expected service life.

### 5.6.1 Manufacturer instructions - electromagnetic emissions

| Emission Test                                             | Compliance |
|-----------------------------------------------------------|------------|
| RF emissions<br>CISPR 11                                  | Group 1    |
| RF emissions<br>CISPR 11                                  | Class B    |
| Harmonic emissions<br>IEC 61000-3-2                       | Class A    |
| Voltage fluctuations / flicker emissions<br>IEC 61000-3-3 | Complies   |

## 5.6.2 Manufacturer instructions - electromagnetic susceptibility

| Immunity Tests                                                                                                  | Compliance level                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                                  | ± 8 kV for contact discharge<br>± 15 kV for air discharge                                                                                                         |
| Radiated RF<br>IEC 61000-4-3<br><br>Proximity fields from RF wireless communications equipment<br>IEC 61000-4-3 | 3 V/m<br>80 MHz – 2,7 GHz<br>80 % AM at 1 kHz<br><br>See Table 1                                                                                                  |
| Fast electrical transients / burst<br>IEC 61000-4-4                                                             | ±2 kV for power line<br>repetition frequency 100 kHz                                                                                                              |
| Surge<br>IEC 61000-4-5                                                                                          | ± 1 kV Line-to-line<br>± 2 kV Line-to-ground                                                                                                                      |
| Conducted RF<br>IEC 61000-4-6                                                                                   | 3 V (0,15 MHz – 80 MHz)<br>6 V in ISM bands between 0,15 MHz and 80 MHz<br>80 % AM at 1 kHz                                                                       |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8                                                      | 30 A/m                                                                                                                                                            |
| Voltage dips, short interruptions on power supply input lines<br>IEC 61000-4-11                                 | 0 % UT; 0,5 cycle<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % UT; 1 cycle and 70 % UT; 25/30 cycle<br>Single phase: at 0°<br>0 % UT; 250/300 cycle |

**Table 1 - IMMUNITY to RF wireless communications equipment**

| Test frequency (MHz)    | Band (MHz)    | Service                                                                  | Modulation                         | Immunity Test Level<br>V/m |
|-------------------------|---------------|--------------------------------------------------------------------------|------------------------------------|----------------------------|
| 385                     | 380 - 390     | TETRA 400                                                                | Pulse modulation 18 Hz             | 27                         |
| 450                     | 430 - 470     | GMRS 460, FRS 460                                                        | FM ± 5 kHz deviation 1<br>kHz sine | 28                         |
| 710<br>745<br>780       | 704 - 787     | LTE band 13, 17                                                          | Pulse modulation 217<br>Hz         | 9                          |
| 810<br>870<br>930       | 800 - 960     | GSM 800/900, TETRA<br>800, iDEN 820, CDMA<br>850,<br>LTE band 5          | Pulse modulation 18 Hz             | 28                         |
| 1 720<br>1 845<br>1 970 | 1 700 - 1 990 | GSM 1800; CDMA<br>1900; GSM 1900;<br>DECT; LTE band 1, 3,<br>4, 25; UMTS | Pulse modulation 217<br>Hz         | 28                         |
| 2 450                   | 2 400 - 2 570 | Bluetooth, WLAN,<br>802.11 b/g/n, RFID<br>2450,<br>LTE band 7            | Pulse modulation 217<br>Hz         | 28                         |
| 5 240<br>5 500<br>5 785 | 5 100 - 5 800 | WLAN 802.11 a/n                                                          | Pulse modulation 217<br>Hz         | 9                          |

**NOTE** There are applied no deviations to requirements of IEC 60601-1-2 ed. 4

**NOTE** There are no known other measures for keeping the basic safety based on EMC phenomena.

**NOTE** Beds equipped with communication module meet standard for IEEE 802.11 b/g/n (2400,0 MHz – 2483,5 MHz, modulation DSSS (IEEE 802.11 b ), OFDM (IEEE 802.11 g/n) 20MHz bandwidth, EIRP = 0,34 W).

## 6 Use and Storage Conditions



### **DANGER!**

#### **Danger to life due to electric shock!**

To ensure the bed's class I protection against electric shocks:

- ▶ Ground the mains.
- ▶ Use exclusively Hospital Grade or Hospital Only receptacles for grounding.

Eleganza 4 is designed for use in rooms for medical purposes. Electrical installations must therefore meet local norms laying down the necessary conditions for electrical installations.

- ▶ Disconnect the bed from the mains in exceptional cases (i.e. lightnings, earthquake).

**Eleganza 4 is not suitable for indoor environments containing flammable gases (except oxygen cylinders).**

## 7 Scope of Delivery and Bed Variants

### 7.1 Delivery

- ▶ Upon receipt, check that the shipment is complete as specified on the delivery note.
- ▶ Notify the carrier and supplier of any deficiencies or damages immediately as well as in writing or make a note on the delivery note.

### 7.2 Scope of Delivery

- Eleganza 4 hospital bed
- Instructions for use

### 7.3 Eleganza 4 Variants

s = standart

o = optional

#### Variable bed features:

- Scales
  - ☐ without scales (without Bed Exit Monitoring) (o)
  - ☐ with scales (with Bed Exit Monitoring or with SafetyMonitor) (s)
- Castors
  - ☐ Tente Integral 150 mm (5.9 in.) single castors (s)
  - ☐ Tente Integral 150 mm (5.9 in.) double castors (o)
  - ☐ Tente Integral 150 mm (5.9 in.) single castors + retractable 5th castor (o)
  - ☐ Tente Integral 150 mm (5.9 in.) double castors + retractable 5th castor (o)
- Control Elements
  - ☐ Nurse Control Panel in both head siderails (o)
  - ☐ iBoard Basic in both head siderails (s)
  - ☐ Attendant Control Panel (s)
  - ☐ Patient Control Panel in both foot siderails (o)
  - ☐ Handset with illuminated buttons and adapter for simple connection Plug and Play (o)
  - ☐ Foot Switch - Bed Height (o)
- Undercarriage Cover
  - ☐ 2-part undercarriage cover (s)
  - ☐ 1-part undercarriage cover (o)
- Linen Shelf (s)
- 1 pair of Mobi-Lift® handles (o)
- Night Light (o)
- Brake Signal (o)
- i-Brake® (o)
  - ☐ without Brake Signal (o)
  - ☐ with Brake Signal (o)
- x-ray cassette holder (o)
- EMR ready bed (o)
- USB (o)
- Safestop (o)
- i-Drive Power (o)
- LINIS SafetyPort
  - ☐ without LINIS SafetyPort (s)
  - ☐ basic hardware preparation for LINIS SafetyPort (CE06: without Integration Module) (o)
  - ☐ complete hardware preparation for LINIS SafetyPort (CE31: with Integration Module) (o)
- Bed ready for Protectors (o)

**NOTE** It is not possible to equip the Eleganza 4 bed with Attendant Control Panel and with Nurse Control Panel at the same time.

## 8 Putting into Service



### **WARNING!**

#### **Risk of injury when working on the bed!**

- ▶ Ensure that the bed is disconnected from the mains connection prior to putting into service, putting out of service and maintenance.
- ▶ Ensure that the castors are locked prior to putting into service, putting out of service and maintenance.



### **CAUTION!**

#### **Material damage due to incorrect putting into service!**

- ▶ Ensure that putting into service is performed exclusively by manufacturer's customer service or trained hospital personnel.



### **CAUTION!**

#### **Material damage due to temperature difference!**

- ▶ If there is a considerable temperature difference between the bed and the place of operation (after transport/storage), leave bed unconnected for 24 hours for the difference to balance itself.

**NOTE** For safe, easy handling, LINET ® recommends having two technicians assemble the bed.

### **Set up the bed as follows:**

- ▶ Unpack the bed.
- ▶ Check the delivery (see Scope of Delivery and Bed Variants).
- ▶ Remove isolating foil from the mains control box (see Accumulator Activation).
- ▶ Install equipment and accessories.
- ▶ In case of delivery with dismantled bed ends, mount the head and foot ends (see Head Board and Foot Board).
- ▶ Set up the bed exclusively on a suitable floor surface (see Transport).
- ▶ Ensure that the mains cable does not collide or get stretched when adjusting the bed.
- ▶ Check that the plug is inserted correctly.
- ▶ Do not leave any extension cords or power strips loose on the floor.
- ▶ Ensure that all of the required mechanical and electrical prevention mechanisms are available on site.
- ▶ There is no mains switch on the bed, i.e. the mains cable is the only means to isolate the bed from the mains. Ensure that the mains cable is always accessible.
- ▶ Have the separable plug of the mains cable changed and maintained exclusively by qualified and trained service technicians authorised by the manufacturer.

## 8.1 Accumulator Activation

### Control Section Placement




Fig. Position of Isolating foil

### Removing the Isolating Foil

To remove isolating foil:

- ▶ Remove isolating foil from mains control box by pulling strap.
- ▶ Check if isolating foil is complete and undamaged as shown.
- ▶ If isolating foil is damaged, contact the manufacturer's service department immediately.

**NOTE** It is recommended to wear gloves when removing the isolating foil.




Fig. Detail of Isolating Foil

## 8.2 Head Board and Foot Board

**Remove head board or foot board as follows:**

- ▶ Unlock Head Board Lock or Foot Board Lock.
- ▶ Pull head board or foot board from sleeve fitting.
- ▶ Lock Head Board Lock or Foot Board Lock.

**Install head board or foot board as follows:**

- ▶ Unlock Head Board Lock or Foot Board Lock.
- ▶ Slide head board or foot board into sleeve fitting.
- ▶ Lock Head Board Lock or Foot Board Lock.




Fig. Head Board Lock

## 8.3 Mattress support platform

Eleganza 4 bed has 4-part Mattress support platform consisting of Backrest, Seat section, Thighrest and Calfrest. The Mattress support platform without x-ray cassette holder has 4 removable covers of Mattress support platform (1, 2, 3 and 4).




Fig. 4-part Mattress support platform

## 8.4 Potential Equalisation

The bed is equipped with a standard protective connector. This connector is used for potential equalisation between the bed and any intravascular or intracardiac device connected to the patient to protect the patient from static electric shocks.

### Use equalisation connector if:


- the patient is connected to any intravascular or intracardiac device.

### Before connecting the patient to an intravascular/intracardiac device:


- ▶ Connect the ground wire of the device to the potential equalisation connector on the bed on which the patient in question is lying.
- ▶ Use a standard hospital connector.
- ▶ Make sure that the connectors match.
- ▶ Make sure that there is no possibility for inadvertent disconnection.

### Before moving the bed:

- ▶ Disconnect the patient from the intravascular or intracardiac device.
- ▶ Disconnect the potential equalisation connector.



*Fig. Potential equalisation connector - female*



*Fig. Potential equalisation - male*

## 8.5 Before Use

### Prepare the bed for use as follows:

- ▶ Connect the bed to the mains.
- ▶ Charge the accumulator.
- ▶ Raise and tilt the mattress support platform to the highest position.
- ▶ Lower and tilt the mattress support platform to the lowest position.
- ▶ Check that the castors as well as main brake work correctly.
- ▶ Check that the bed extension works correctly.
- ▶ Check that it is possible to remove the head and foot boards.
- ▶ Check all of the functions on the control elements.
- ▶ Check that the siderails function properly.
- ▶ Dispose of all packaging (see Disposal).

## 8.6 Transport

### For a safe transport, observe the following:

- ▶ Ensure that no cables are run over when moving a bed.
- ▶ Ensure that the mains cable is attached with a hook (at the head end of the bed).
- ▶ Ensure that the castors are unlocked before moving the bed during the loading/unloading process (see Castor Control).
- ▶ Adjust bed height to at least 20 cm below maximum height.
- ▶ Push bed by handles on head or foot end.
- ▶ Move the bed exclusively on suitable floor surfaces.
- ▶ For longer distances, ensure that the castor steering function is activated.
- ▶ Ensure that the brakes are released while moving the bed.

### Suitable surfaces:

- Tile
- Hard linoleum
- Poured flooring

### Unsuitable surfaces:

- Too soft, unsealed or defective flooring
- Soft wooden flooring
- Soft and porous stone floors
- Carpeted floors with underlay
- Soft linoleum

## 8.7 Firmware

The bed includes firmware that can be updated only by an authorised service technician.

This firmware is protected against unauthorised access by mechanical housing (tool is needed to access), by seal (components with processor are sealed), by exclusive compatibility with an authorised software tool and by check of compatibility of the new firmware with the bed.

## 9 Power Cable

Attachment plug is means of connecting and disconnecting bed from the mains.  
Mains power cable must be attached with a hook at the head end of the bed during transport.



### CAUTION!

**Disconnecting bed from the mains does not stop motions of the bed!**

- ▶ Stop the bed before disconnection bed from the mains.

**Where the integrity of the external protective conductor in the installation or its arrangement is in doubt**

- ▶ operate the bed from internal accumulator only.

## 10 Accumulator



### WARNING!

**When the bed is not connected to the mains and accumulator is not sufficiently charged all electrical functions of the bed are blocked!**

### Purpose

**The accumulator serves as a backup during power failures, during transport of patient or for emergency bed positioning.**

- ▶ Use only accumulator approved by the manufacturer.
- ▶ Check the functionality of the accumulator at least once a month in accordance with the instructions for use and have the accumulator changed if necessary.

**The manufacturer will assume no responsibility for any damage to the bed or the accumulator caused by:**

- non-observance of the manufacturer's instructions in the instructions for use
- using accumulator not approved by the manufacturer

### Warranty

The manufacturer provides a **6-month** warranty for the full function of the accumulator.

**Accumulator lifetime could be up to 5 years if operated under optimum conditions.**

**Accumulator capacity can be significantly reduced if:**

- ▶ too high ambient temperature
- ▶ many accumulator charge/discharge cycles
- ▶ recurrence of deep discharge
- ▶ bed is often powered only by the accumulator

### Charging

The accumulator supplied with the bed is delivered insufficiently charged.

Charging of the accumulator before use of the bed takes approximately 4 hours.

Charging of the accumulator runs automatically when the power cable is connected to the mains.

**To charge the accumulator:**

- ▶ Connect the bed to the mains.

### Storage

**For declared lifetime period of leaded accumulator is recommended during storage:**

- ▶ To prevent accumulator from deep discharging and to keep accumulator at least partly charged by regular recharging
- ▶ To store accumulator on the places with temperature from 10°C to 40°C
- ▶ To prevent accumulator from being in the sunshine

## Signalisation

The LED (on Nurse Control Panel or Attendant Control Panel) indicates the accumulator charge status:

| Yellow LED                                                    | Accumulator charge status                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Not lit                                                       | Accumulator capacity is sufficient (charging completed)                                                                                                                                                                                                                                                                                     |
| Short flashing (shortly lit, longer not lit) (circa 1.8 sec.) | Accumulator is charging - continue charging until the LED is extinguished. In emergency cases, the accumulator can be used as a backup power source for a short period. If LED is still flashing after 12 hours of charging or stops flashing, but you can not position with bed, accumulator is defective or broken. Contact manufacturer. |
| Short flashing (0,2s lit, 0,2s not lit)                       | Only CPR function can be used.                                                                                                                                                                                                                                                                                                              |
| Long flashing (longer lit, shortly not lit) (circa 0.2 sec.)  | Low accumulator voltage - accumulator can not be used as a backup power supply even for a short period; accumulator is completely discharged or defective (if this type of signalisation persists, it is necessary to replace the accumulator - service action)                                                                             |
| Lit continuously, when bed is connected to the mains.         | Accumulator absence or failure condition (accumulator is connected incorrectly, line between the power supply and accumulator is broken or accumulator fuses are faulty); contact service department of the manufacturer in case of such signalisation.                                                                                     |

## 10.1 Replacing the accumulator



### CAUTION!

#### Damage to the bed due to incorrect accumulator replacement!

- Have the accumulator replaced exclusively by qualified personnel.
- Exclusively use accumulator approved by the manufacturer.



### CAUTION!

#### Material damage due to overheating!

If the accumulator is faulty, degassing may occur. In rare cases this might cause deformations of the accumulator case, control panel housing or cable.

- Stop using the bed immediately (see Removing the Bed from Service).
- Inform the manufacturer's service department.



### CAUTION!

#### Risk of reducing accumulator durability due to incorrect use!

- Use bed on accumulator only in crisis situations (e.g.: power blackout, patient complications during transport, etc.).
- After reconnecting bed to the mains charge accumulator to full capacity (see chart Accumulator charge status).
- Have accumulator replaced exclusively by a qualified service organisation.
- The manufacturer recommends to replace the accumulator by qualified service organization after 2 years of use.

After this period the supposed service life of accumulator ends and the manufacturer cannot guarantee the accumulator service life after this period.

### Faulty accumulator

This status is indicated by the accumulator status indicator being constantly lit.

The accumulator is regarded as faulty if at least one of the following conditions applies:

- Accumulator charging constantly
- Low voltage on accumulator
- Low charging current of accumulator

**NOTE** These statuses are summarised and written to black box in control unit.

### To cancel this status:

- Press button .

## Discharged accumulator

This status will be cancelled automatically when the bed switches to sleep mode.

The accumulator is regarded as discharged if the following condition is met:

- ▶ Defined decrease of voltage depending on discharging current
- ▶ This status is indicated by the accumulator status indicator flashing quickly.
- ▶ It is possible to use the electric CPR only to position the bed.

To cancel this status:

- ▶ Press button  .

## 10.2 Removing the Bed from Use

Remove the bed from use as follows:

- ▶ Disconnect the bed from the mains.
- ▶ Disconnect the ground wire.
- ▶ Deactivate the accumulator.
- ▶ Remove accessories.

To prevent damage during storage:

- ▶ Pack or cover the bed and accessories.
- ▶ Ensure that storage conditions are the same as the operating conditions.

## 10.3 Deactivating the Accumulator

To avoid damaging the bed and the environment during storage:

- ▶ Deactivate the accumulator.

To deactivate the accumulator:

- ▶ Disconnect the bed from the mains.
- ▶ Disconnect the ground wire.
- ▶ Activate the keypad by pressing the button  .
- ▶ Press the Thighest Up + Thighest Down + Trendelenburg Tilt buttons at the same time and hold them for 3 seconds. The accumulator is deactivated.




Fig. Deactivation of the accumulator (Nurse Control Panel)




Fig. Deactivation of the accumulator (Attendant Control Panel)

To activate the accumulator again:

- ▶ Connect Power Cable to the mains.

## 11 Manipulation



### **WARNING!**

#### **Risk of injury when adjusting the bed!**

- ▶ Ensure that there are no body parts between the mattress support platform elements and the mattress support platform frame when adjusting the bed.
- ▶ Ensure that there are no body parts below the mattress support platform frame before adjusting the bed.

### 11.1 Siderails

The split siderails are components of the bed in contact with patient. A pneumatic spring supports the operation of the split siderails. The nursing personnel are responsible for the siderails being raised up while the patient is in bed.



### **WARNING!**

#### **Risk of injury, damaging or unintentional movement of the bed due to incorrect placement of accessories or Handset!**

- ▶ Never place any accessories or Handset on the siderails in the area where keyboards are located.
- ▶ Never place Handset on the edge of siderail.

The correct placement of Handset is shown at following pictures.



### **WARNING!**

#### **Risk of injury due to incorrectly latched siderail!**

- ▶ Ensure that siderail is secured in the upper or lower position.



### **WARNING!**

#### **Risk of injury due to incorrect position of siderails!**

- ▶ Ensure that siderails are folded up while the patient is in bed.

## SIDERAIL DESCRIPTION (version with Nurse Control Panel and iBoard Basic)




Fig. Siderails (version with Nurse Control Panel and iBoard Basic)

- 1. Nurse Control Panel
- 2. iBoard Basic
- 3. Siderail Handle
- 4. Correct Placement of Handset
- 5. Siderail Release Handle
- 6. Angle Indicator

### MANIPULATION

#### To raise siderails up:

- ▶ Grab siderail by Siderail Handle (3).
- ▶ Pull siderail up until it latches. You will hear audible „click“.

#### To release siderails down:

- ▶ Grab siderail by Siderail Handle (3).
- ▶ Unlock siderail by pulling Siderail Release Handle (5) to yourself.
- ▶ Fold down siderail slowly.

## SIDERAIL DESCRIPTION (version with Nurse Control Panel)




Fig. Siderails (version with Nurse Control Panel)

- 1. Nurse Control Panel
- 2. Siderail Handle
- 3. Correct Placement of Handset
- 4. Siderail Release Handle
- 5. Angle Indicator

### MANIPULATION

#### To raise siderails up:

- ▶ Grab siderail by Siderail Handle (2).
- ▶ Pull siderail up until it latches. You will hear audible „click“.

#### To release siderails down:

- ▶ Grab siderail by Siderail Handle (2).
- ▶ Unlock siderail by pulling Siderail Release Handle (4) to yourself.
- ▶ Fold down siderail slowly.

## SIDERAIL DESCRIPTION (version with iBoard Basic)




Fig. Siderails (version with iBoard Basic)

- 1. iBoard Basic
- 2. Siderail Handle
- 3. Correct Placement of Handset
- 4. Siderail Release Handle
- 5. Angle Indicator

## MANIPULATION

### To raise siderails up:

- ▶ Grab siderail by Siderail Handle (2).
- ▶ Pull siderail up until it latches. You will hear audible „click“.

### To release siderails down:

- ▶ Grab siderail by Siderail Handle (2).
- ▶ Unlock siderail by pulling Siderail Release Handle (4) to yourself.
- ▶ Fold down siderail slowly.

## SIDERAIL DESCRIPTION (version without control panels)




Fig. Siderails (version without control panels)

- 1. Siderail Handle
- 2. Correct Placement of Handset
- 3. Siderail Release Handle
- 4. Angle Indicator

### MANIPULATION

#### To raise siderails up:

- ▶ Grab siderail by Siderail Handle (1).
- ▶ Pull siderail up until it latches. You will hear audible „click“.

#### To release siderails down:

- ▶ Grab siderail by Siderail Handle (1).
- ▶ Unlock siderail by pulling Siderail Release Handle (3) to yourself.
- ▶ Fold down siderail slowly.

## 11.2 Castor Control



### CAUTION!

#### Material damage due to incorrect transport and involuntary movement!

- ▶ Prior to transport, ensure that the bed is disconnected from the mains.
- ▶ Ensure that the castors are braked prior to assembly, disassembly and maintenance.
- ▶ Ensure that the castors are braked when the bed is occupied.
- ▶ Hang the mains cable on the transport hook on the bed during transport.
- ▶ Have the bed transported exclusively by nursing personnel and by at least 2 persons.



### CAUTION!

#### Minimal clearance underneath the bed (standard version with 15 cm castors) is 11,3 cm!

- ▶ Observe the path for any obstacles and avoid collisions and possible damages of any bed's part on the undercarriage.
- ▶ Do not use bed lifts and hoists for lifting the bed.

The bed is equipped with central castor's control and brake system.  
The control levers are located in the four corners of the undercarriage.




Fig. Positions of Castor Control Lever

### Castor control lever positions:

#### 1. Forward Movement - Steering (GREEN PEDAL DOWN)

An arrested castor determines the direction of movement. The bed moves straight ahead. If the bed is equipped with a fifth castor, this castor determines the direction of movement.

#### 2. Unrestricted Movement

All four castors are unlocked.

#### 3. Braked (RED PEDAL DOWN)

All four castors are braked.

## 11.3 CPR Backrest Release



### WARNING!

**Risk of injury due to lowering the backrest too quickly!**

- ▶ Ensure that the siderails are in the low position.
- ▶ Ensure that there are no body parts between any movable parts of the bed.
- ▶ Push the Backrest down using the mattress guard handle only.

The bed allows quick, mechanical lowering of the backrest for emergency procedures (CPR).

### Set the CPR position as follows:

- ▶ Pull and hold release handle.
- ▶ Push Backrest down.




Fig. CPR lever

## 11.4 Control Elements

The bed is operated by different control elements.

### Control elements depending on the model:

- Nurse Control Panel in both head siderails
- iBoard Basic in both head siderails
- Attendant Control Panel
- Handset with illuminated buttons and with adapter for easy connection (Plug and Play)
- Patient Control Panel in both foot siderails
- Foot Switch - Bed Height

**Disabling individual functions on the Attendant Control Panel or on the Nurse Control Panel will affect all control elements.**

### If the bed does not react to individual position settings:

- Check whether the function is disabled on the Attendant Control Panel or on the Nurse Control Panel.

| POSITIONING                               | Nurse Control Panel | Attendant Control Panel | Handset | Patient Control Panel | Foot Switch- Bed Height |
|-------------------------------------------|---------------------|-------------------------|---------|-----------------------|-------------------------|
| Backrest                                  | ✓                   | ✓                       | ✓       | ✓                     |                         |
| Thighrest                                 | ✓                   | ✓                       | ✓       | ✓                     |                         |
| Bed Height                                | ✓                   | ✓                       | ✓       |                       | ✓                       |
| Autocontour                               | ✓                   | ✓                       | ✓       | ✓                     |                         |
| Examination Position                      |                     | ✓                       |         |                       | ✓                       |
| Emergency Trendelenburg Position          | ✓                   | ✓                       |         |                       |                         |
| Anti-Trendelenburg and Trendelenburg Tilt | ✓                   | ✓                       |         |                       |                         |
| CPR Position                              | ✓                   | ✓                       |         |                       |                         |
| Cardiac Chair Position                    | ✓                   | ✓                       |         |                       |                         |
| Mobilization Position                     | ✓                   |                         |         |                       |                         |

## 11.4.1 Nurse Control Panel (optional)

The Nurse Control Panel is the main Control Element for the caregivers. It is integrated in the outside of both head siderails.

- Ensure that exclusively trained nursing staff operates the Nurse Control Panel.




Fig. Nurse Control Panel

1. Mobilization Position Button
2. Trendelenburg Tilt Button
3. Anti-Trendelenburg Tilt Button
4. LIGHT CONTROL Button
5. Autocontour Adjustment Buttons (simultaneous movement of the Backrest and Thighrest)
6. LOCK Button
- 6a. Locked Backrest LED
- 6b. Locked Thighrest LED
- 6c. Locked Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt LED
7. Central STOP Button
8. Backrest Adjustment Buttons
9. Thighrest Adjustment Buttons
10. Mains Power LED
11. Accumulator LED
12. Bed Height Adjustment Buttons
13. GO Button
14. Emergency Trendelenburg Position Button
15. CPR (Resuscitation) Position Button
16. Cardiac Chair Position Button

Positioning buttons 1, 2, 3, 5, 8, 9, 12, 14, 15 and 16 are explained in chapter **Bed Positioning**.

### GO BUTTON

The button  activates the keyboard of all Control Elements.

Pressing button  will keep the keyboard active for 3 minutes.

Pressing a button will keep the keyboard active for another 3 minutes.

### During this time the following is possible:

- Adjusting individual Mattress support platform elements by pressing the corresponding positioning buttons.
- Disabling individual functions with the lock buttons.

**NOTE:** To activate CPR function (button ) the button  is not needed.

## STOP BUTTON

Pressing button  immediately stops all electronic bed movements.

## LOCK

**To lock Backrest Adjustment:**


- ▶ Press button .

LED Lock **6a**, **6b** and **6c** are flashing.

- ▶ During flashing of the LEDs press any button **8**.

LED Lock **6a** is lit.

Backrest Adjustment is disabled using any Control Element.



**To lock Thighrest:**


- ▶ Press button .

LED Lock **6a**, **6b** et **6c** are flashing.

- ▶ During flashing of the LEDs press any button **9**.

LED Lock **6b** is lit.

Thighrest is disabled using any Control Element.



**To lock Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt :**


- ▶ Press button .

LED Lock **6a**, **6b** et **6c** are flashing.

- ▶ During flashing of the LEDs press any button **12**.

LED Lock **6c** is lit.

Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt are disabled using any Control Element.



## UNLOCK

**To unlock Backrest Adjustment:**


- ▶ Press button .

LED Lock **6a**, **6b** et **6c** are flashing.

- ▶ During flashing of the LEDs press any button **8**.

LED Lock **6a** is unlit.

Backrest Adjustment is enabled again.



**To unlock Thighrest:**


- ▶ Press button .

LED Lock **6a**, **6b** et **6c** are flashing.

- ▶ During flashing of the LEDs press any button **9**.

LED Lock **6b** is unlit.

Thighrest is enabled again.



To unlock **Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt:**


- Press button .

LED Lock **6a, 6b** et **6c** are flashing.

- During flashing of the LEDs press any button **12**.

LED Lock **6c** is unlit.

Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt are enabled again.



**LOCKED FUNCTION SIGNALISATION**

If **LED 6a** is lit, Backrest Adjustment is locked.

If **LED 6a** is not lit, Backrest Adjustment is unlocked.

If **LED 6b** is lit, Thighrest is locked.

If **LED 6b** is not lit, Thighrest is unlocked.

If **LED 6c** is lit, Bed Height and Trendelenburg Tilt and Anti-Trendelenburg Tilt are locked.

If **LED 6c** is not lit, Bed Height and Trendelenburg Tilt and Anti-Trendelenburg Tilt are unlocked.

**ACCUMULATOR INDICATOR**

Signalisation of Accumulator LED  is described in chapter **Accumulator**.

**MAINS POWER LED** 

| Status       | Meaning                     |
|--------------|-----------------------------|
| lit LED      | connected to the mains      |
| unlit LED    | disconnected from the mains |
| flashing LED | system error                |

**LIGHT CONTROL BUTTON**

To turn off lights of the control elements:

- Press button .

To turn on lights of the control elements again:

- Press button .

## 11.4.2 iBoard Basic (standard)

The iBoard Basic is the standard control element for the caregivers. It is integrated in the outside of both head siderails. Only bed version with scales can be equipped with iBoard Basic. The iBoard Basic serves for control of scales and control of Bed Exit Monitoring.

- Ensure that exclusively trained nursing staff operates the iBoard Basic.




Fig. iBoard Basic

### Scales Section (only version with scales)

Eleganza 4 is optionally equipped with a weighing system that allows weighing the patient in bed. There are control buttons and display for the weighing system on the Scales section of iBoard Basic. Scales functions are described in chapter **Scales Control**.

### Bed Exit Monitoring Section (only version with scales)

Eleganza 4 is optionally equipped with a Bed Exit Monitoring system that monitors patient's presence in bed and triggers alarms when patient is not present in bed. There are control buttons and display for the Bed Exit Monitoring on the SafetyMonitor section of iBoard Basic. Functions of Bed Exit Monitoring are described in chapter **Bed Exit Monitoring**.

## Statuses (iBoard Basic)

| Signalisation                                                                                                                                                                                                   | Meaning                                                                                                              | Required Action                                                                                                        |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                | Function locked.                                                                                                     | Unlock function if required!                                                                                           |
|                                                                                                                                | GO Button not activated.                                                                                             | Press GO Button!                                                                                                       |
| 0°                                                                                                                                                                                                              | Optionally: Horizontal position was reached during tilting.                                                          | Press corresponding button to continue in positioning.                                                                 |
| <b>OVERLOAD</b> +                                                                                                              | Safe Working Load exceeded (more than 10 kg over Safe Working Load).                                                 | Remove load!                                                                                                           |
| <b>STOP SERVICE</b> +  +                      | System Fatal Error.                                                                                                  | Contact service department approved by manufacturer.                                                                   |
| <b>SCALE / BEA DISCONNECTED</b> +  +         | Scale module disconnected and Bed Exit monitoring disabled.                                                          | Contact service department approved by manufacturer.                                                                   |
| <b>SAFE STOP</b> +                                                                                                           | Movement of the Mattress support platform stopped by function Safestop.                                              | Remove an object from undercarriage to continue in adjusting the bed height.                                           |
|                                                                                                                              | Insufficient load for Bed Exit monitoring.                                                                           | Place patient on the bed to enable the Bed Exit monitoring.                                                            |
|                                                                                                                              | Disconnected from the mains power supply.                                                                            | Connect bed to the mains power.                                                                                        |
| <b>BED EXIT</b> + <br><br><b>ALARM</b> +  | The patient has left the bed (Outer Zone monitoring) or the patient has left the Inner Zone (Inner Zone monitoring). | Check the patient and turn off the Bed Exit Alarm.                                                                     |
| <b>HIGH</b>                                                                                                                                                                                                     | Safe Working Load exceeded (from 4,5 kg to 10 kg over Safe Working Load).                                            | Remove load!                                                                                                           |
| <b>LOW</b>                                                                                                                                                                                                      | The bed is underloaded.                                                                                              | Ensure the mattress support platform is not lifted inconveniently by something and the scales are appropriately tared. |

### 11.4.3 Attendant Control Panel

The Attendant Control Panel is a standard control element. The Attendant Control Panel can be hung on the foot board or on the siderails if required. The Attendant Control Panel can be stored in the Linen Shelf. It is possible to hold the Attendant Control Panel in the hand while operating.

- Ensure that exclusively trained nursing staff operates the Attendant Control Panel.




Fig. Attendant Control Panel

1. GO Button
2. Thighrest Adjustment Button
3. Thighrest Adjustment Lock Button and LED
4. Backrest Adjustment Button
5. Backrest Lock Button and LED
6. Bed Height Adjustment Button
7. Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt Lock Button and LED
8. Autocontour Adjustment Button
9. Foot Switch Lock Button
10. Foot Switch Lock LED
11. Trendelenburg Tilt Button
12. Anti-Trendelenburg Tilt Button
13. Examination Position Button
14. CPR Position Button
15. Emergency Trendelenburg Position Button
16. Cardiac Chair Position Button
17. Accumulator Charge Status LED
18. Mains Power LED
19. LIGHT CONTROL Button
20. STOP Button

#### Central STOP Button



The central STOP Button immediately interrupts all bed movements in case of unauthorized bed positioning or an electronic failure. Pressing the central STOP Button for at least 0.3 seconds immediately stops all electronic bed movements.

## Activating GO Button



The GO Button activates the keypad of all control elements for 3 minutes. A GO Button is included on a number of different control elements. The function of the GO Button is identical on all control elements.

### During this time the following is possible:

- ▶ Adjusting individual mattress support platform elements by pressing the corresponding function buttons.
- ▶ Disabling individual functions with the lock buttons.

Each time a function button is pressed, the keypad will remain active for another 3 minutes.

## Function Buttons

The function buttons **2, 4, 6, 8, 11, 12, 13, 14, 15** and **16** are described in chapter **Bed Positioning**.

**NOTE** Pressing two function buttons at the same time will be recognized as an error by the controller. The controller will interrupt immediately all bed movements and display shows adequate pop up.

## Lock

### To lock Backrest Adjustment:

- ▶ Press button **5**.

Corresponding LED on Attendant Control Panel and on Nurse Control Panel (**6a**) are lit.

Backrest Adjustment is disabled using any Control Element.

### To lock Thighrest Adjustment:

- ▶ Press button **3**.

Corresponding LED on Attendant Control Panel and on Nurse Control Panel (**6b**) are lit.

Thighrest Adjustment is disabled using any Control Element.

### To lock Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt Adjustment:

- ▶ Press button **7**.

Corresponding LED on Attendant Control Panel and on Nurse Control Panel (**6c**) are lit.

Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt Adjustment are disabled using any Control Element.

### To lock Foot Switches:

- ▶ Press button **9**.

Corresponding LED (**10**) on Attendant Control Panel is lit.

Bed Height and Examination Position are disabled using Foot Switches.

## Unlock

### To unlock Backrest Adjustment:

- ▶ Press button **5**.

Corresponding LED on Attendant Control Panel and on Nurse Control Panel (**6a**) are not lit.

Backrest Adjustment is enabled again.

### To unlock Thighrest Adjustment:

- ▶ Press button **3**.

Corresponding LED on Attendant Control Panel and on Nurse Control Panel (**6b**) are not lit.

Thighrest Adjustment is enabled again.

### To unlock Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt Adjustment:

- ▶ Press button **7**.

Corresponding LED on Attendant Control Panel and on Nurse Control Panel (**6c**) are not lit.

Bed Height, Trendelenburg Tilt and Anti-Trendelenburg Tilt Adjustment are enabled again.

### To unlock Foot Switches:

- ▶ Press button **9**.

Corresponding LED (**10**) on Attendant Control Panel is no lit.

Bed Height and Examination Position are enabled again.

## MAINS POWER LED



| Status       | Meaning                     |
|--------------|-----------------------------|
| lit LED      | connected to the mains      |
| unlit LED    | disconnected from the mains |
| flashing LED | system error                |

### 11.4.4 Handset (optional)

The handset is available with illuminated keyboard.

**NOTE** The illumination is activated for 7s if any button was pressed and the illumination is activated for 3 minutes if GO Button was pressed.




Fig. Handset

The function buttons **1**, **3**, **5** and **8** are described in chapter **Bed Positioning**.

#### To switch on the flashlight:

- Press flashlight button **6**.

**NOTE** Depending on the patient's condition, the nursing staff decides whether the patient is allowed to adjust the bed's position.

#### If required, prevent the patient from adjusting the bed as follows:

- Disable functions.

**NOTE** An adapter for the handset is available. The adapter enables quick installation and removal (e.g. replacing a defective handset, using the handset for another bed).

## 11.4.5 Foot Switch - Bed Height (optional)

The foot control is optional and allows setting the Height of the bed or Examination Position with one's feet. Press the selected pedal twice in 3 seconds. Bed Height Foot Control is activated for 20s after this procedure.





Fig. Foot Switch - Bed Height

The use of Foot Switch - Bed Height is described in the chapter **Bed Positioning**.

## 11.4.6 Patient Control Panels

The patient control panels integrated in the foot siderails allow the patient to adjust the positions of the Backrest, Thighrest and Autocontour.




**NOTE** Functions on the Patient Control Panel in the foot siderails are disabled when the foot siderail is in lower position.

Fig. Patient Control Panel in the foot siderails

# 11.5 Bed Positioning

## 11.5.1 Backrest



To position Backrest use:


- ▶ Nurse Control Panel
- ▶ Attendant Control Panel
- ▶ Handset
- ▶ Patient Control Panel



*Fig. Backrest Angle on iBoard Basic Display*

iBoard Basic Display shows Backrest Angle.

During continuous positioning Backrest stops automatically in 30 and 45 degrees. To continue in positioning release the button, then press and hold button until the desired position is reached.



1. Backrest Up

2. Backrest Down

*Fig. Backrest Adjustment Button  
(Nurse Control Panel, Patient  
Control Panel)*

**Nurse Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Backrest Adjustment Button until intended position is reached.

**Patient Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Backrest Adjustment Button until intended position is reached.

1. Backrest Up

2. Backrest Down




Fig. Backrest Adjustment Button  
(Attendant Control Panel)

**Attendant Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Backrest Adjustment Button until intended position is reached.





Fig. Backrest Adjustment Button  
(Handset)

**Handset:**


- ▶ Press button .
- ▶ Press selected part of Backrest Adjustment Button until intended position is reached.

## 11.5.2 Thighrest



**To position Thighrest use:**

- ▶ Nurse Control Panel
- ▶ Attendant Control Panel
- ▶ Handset
- ▶ Patient Control Panel



*Fig. Thighrest Adjustment Button (Nurse Control Panel, Patient Control Panel)*

**Nurse Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Thighrest Adjustment Button until intended position is reached.

**Patient Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Thighrest Adjustment Button until intended position is reached.

**1. Thighrest Up**


**2. Thighrest Down**



*Fig. Thighrest Adjustment Button (Attendant Control Panel)*

**Attendant Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Thighrest Adjustment Button until intended position is reached.




*Fig. Thighrest Adjustment Button (Handset)*

**Handset:**

- ▶ Press button .
- ▶ Press selected part of Thighrest Adjustment Button until intended position is reached.

## 11.5.3 Mechanical Calfrest Adjustment



### **WARNING!**

**Incorrect manipulation can injure the patient or the user!**


- ▶ Hold the Calfrest by both Calfrest Handles at all times when lowering the Calfrest.
- ▶ Lower the Calfrest carefully to prevent it from falling suddenly.

### **To raise the Calfrest:**

- ▶ lift Calfrest by the handle to intended position.
- ▶ lower the Calfrest so that catch fits in the ratchet-bar.


### **To lower the Calfrest:**

- ▶ lift Calfrest slightly by the handle.
- ▶ lower Calfrest to intended position.
- ▶ ensure the catch fits in the ratchet-bar.



*Fig. Position of Calfrest Handle*

## 11.5.4 Bed Height



To position Bed Height use:

- ▶ Nurse Control Panel
- ▶ Attendant Control Panel
- ▶ Handset
- ▶ Foot Switch - Bed Height
- ▶ Mobi-Lift Handle

The bed stops in the low position during lowering, mattress support platform is adjusted to flat position and beep is performed. After this adjustment the lowering continues to the lowest position.




Fig. Bed Height Adjustment Button (Nurse Control Panel)

Nurse Control Panel:

- ▶ Press button .
- ▶ Press selected part of Bed Height Adjustment Button until intended position is reached.

1. Mattress support platform Up
2. Mattress support platform Down




Fig. Bed Height Adjustment Button  
(Attendant Control Panel)

#### Attendant Control Panel:

- ▶ Press button .
- ▶ Press selected part of Bed Height Adjustment Button until intended position is reached.




Fig. Bed Height Adjustment Button  
(Handset)

#### Handset:

- ▶ Press button .
- ▶ Press selected part of Bed Height Adjustment Button until intended position is reached.

1. Mattress support platform Up
2. Mattress support platform Down




Fig. Bed Height Adjustment Pedals (Foot Switch - Bed Height)

#### Foot Switch - Bed Height:

- ▶ Press and release the selected Bed Height Adjustment Pedal to activate the panel.
- ▶ Press and hold selected Bed Height Adjustment Pedal until intended position is reached.

**NOTE** Press the selected pedal twice in 3 seconds. Bed Height Foot Control is activated for 20s after this procedure.





Fig. Mobi-Lift Bed Height Button

#### Mobi-Lift:

- ▶ Press button .
- ▶ Press selected part of Mobi-Lift Bed Height Button until intended position is reached.

## 11.5.5 Autocontour



To position Autocontour use:

- ▶ Nurse Control Panel
- ▶ Attendant Control Panel
- ▶ Handset
- ▶ Patient Control Panel




Fig. Autocontour Adjustment Button  
(Nurse Control Panel, Patient Control Panel)

**Nurse Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Autocontour Adjustment Button until intended position is reached.

**Patient Control Panel:**

- ▶ Press button .
- ▶ Press selected part of Autocontour Adjustment Button until intended position is reached.

**1. Autocontour Up**

**2. Autocontour Down**




Fig. Autocontour Adjustment Button (Attendant Control Panel)

## Attendant Control Panel:

- ▶ Press button .
- ▶ Press selected part of Autocontour Adjustment Button until intended position is reached.

### 1. Autocontour Up

### 2. Autocontour Down





Fig. Autocontour Adjustment Button (Handset)

## Handset:

- ▶ Press button .
- ▶ Press selected part of Autocontour Adjustment Button until intended position is reached.

## 11.5.6 Examination Position




### To position Examination Position use:

- ▶ Attendant Control Panel
- ▶ Foot Switch - Bed Height



*Fig. Examination Position Button (Attendant Control Panel)*



*Fig. Examination Position Pedal (Foot Switch - Bed Height)*

## Attendant Control Panel:

- ▶ Press button  .
- ▶ Press Examination Position Button until intended position is reached.

## Foot Switch - Bed Height:

- ▶ Press and release the middle pedal to activate the panel.
- ▶ Press and hold Examination Position Pedal until intended position is reached.

**NOTE** There must be less than 20 seconds between these 2 steps.

## 11.5.7 Emergency Trendelenburg Position




Fig. Emergency Trendelenburg Position/Trendelenburg Tilt

To position Emergency Trendelenburg Position use:

- ▶ Nurse Control Panel
- ▶ Attendant Control Panel

**Trendelenburg position provides anti-shock conditions for the patient.**  
**During Trendelenburg Position Mattress support platform is straightened in the tilt.**  
**The LOCK function does not influence the Emergency Trendelenburg Position function!**




Fig. Emergency Trendelenburg Position Button (Nurse Control Panel)




Fig. Emergency Trendelenburg Position Button (Attendant Control Panel)

**Nurse Control Panel:**

- ▶ Press button .
- ▶ Press Emergency Trendelenburg Position Button until intended position is reached.

**Attendant Control Panel:**

- ▶ Press button .
- ▶ Press Emergency Trendelenburg Position Button until intended position is reached.

## 11.5.8 Anti-Trendelenburg Tilt and Trendelenburg Tilt




Fig. Anti-Trendelenburg Tilt

To position Trendelenburg Tilt or Anti-Trendelenburg Tilt use:

- ▶ Nurse Control Panel
- ▶ Attendant Control Panel

Depending on configuration, the bed stops in the horizontal position (0°) during tilting.  
To continue tilting press the corresponding tilt button.




Fig. Trendelenburg Tilt Button  
(Nurse Control Panel)

**Nurse Control Panel:**

- ▶ Press button .
- ▶ Press Trendelenburg Tilt Button until intended position is reached.





Fig. Anti-Trendelenburg Tilt Button  
(Nurse Control Panel)

**Nurse Control Panel:**

- ▶ Press button .
- ▶ Press Anti-Trendelenburg Tilt Button until intended position is reached.



*Fig. Trendelenburg Tilt Button  
(Attendant Control Panel)*

**Attendant Control Panel:**

- ▶ Press button .
- ▶ Press Trendelenburg Tilt Button until intended position is reached.




*Fig. Anti-Trendelenburg Tilt  
Button (Attendant Control  
Panel)*

**Attendant Control Panel:**

- ▶ Press button .
- ▶ Press Anti-Trendelenburg Tilt Button until intended position is reached.

## 11.5.9 CPR Position



To position CPR Position use:

- Nurse Control Panel
- Attendant Control Panel

In CPR Position the mattress support platform is adjusted to the low position and all the parts of the mattress support platform are in the basic (flat) position.



*Fig. CPR Position Button (Nurse Control Panel)*

**Nurse Control Panel:**

- Press CPR Position Button until intended position is reached.



*Fig. CPR Position Button (Attendant Control Panel)*

**Attendant Control Panel:**

- Press CPR Position Button until intended position is reached.

## 11.5.10 Mechanical Bed Extension





Fig. Position of Bed Extension Safety Catch

**There are 2 positions of the Bed Extension.**

### **To extend the bed:**

- ▶ Pull out the Safety Catches at Foot End on the both sides of the frame.
- ▶ Turn the Safety Catches by 90°.
- The Safety Catches are released.
- ▶ Pull the Foot Board to the selected position.
- ▶ Turn the Safety Catches by 90° and push them in to fix the extension position.

### **To shorten the bed:**

- ▶ Pull out the Safety Catches at Foot End on both sides of the frame.
- ▶ Turn the Safety Catches by 90°.
- The Safety Catches are released.
- ▶ Push the Foot Board to the selected position.
- ▶ Turn the Safety Catches by 90° and push them in to fix the extension position.

## Mattress Holder

It is possible to adjust Mattress Holder when the bed without mattress is extended.

To adjust the Mattress Holder:

- ▶ Remove the Calfrest cover from the mattress support platform.
- ▶ Insert the catches on both sides of the Mattress Holder to the selected positions.
- ▶ Put the Calfrest cover back to its place on the mattress support platform.







Fig. Adjustment of the Mattress Holder



*Fig. Basic position of the Mattress Holder*




*Fig. First extension of the Mattress Holder*



*Fig. Second extension of the Mattress Holder*

## 11.5.11 Cardiac Chair Position



**To position Cardiac Chair Position use:**

- Nurse Control Panel
- Attendant Control Panel



*Fig. Cardiac Chair Position Button (Nurse Control Panel)*

**Nurse Control Panel:**

- Press button .
- Press Cardiac Chair Position Button until intended position is reached.




*Fig. Cardiac Chair Position Button (Attendant Control Panel)*

**Attendant Control Panel:**

- Press button .
- Press Cardiac Chair Position Button until intended position is reached.

## 11.5.12 Mobilisation Position



To position Mobilisation Position use:

- Nurse Control Panel

In Mobilisation Position bed is descending to the lowest Bed Height and Backrest reaches the maximum angle.



*Fig. Mobilisation Position  
Button (Nurse Control Panel)*

**Nurse Control Panel:**

- Press button .
- Press Mobilisation Position Button until intended position is reached.

## 11.5.13 Ergoframe

Ergoframe® is the kinematic system of Backrest and Thighrest Adjustment resulting in extension of the Mattress support platform in the seat section.

Ergoframe® enlarges the space for pelvic area during Auto-contour. Because of increasement of the space the force applied results in decrease of the pressure that can cause pressure injuries in the pelvic area.

Ergoframe maintains a stable ergonomic position of the body and spine of the patient, thus limiting unwanted movement of the patient by moving down or up in beds. Unified movement eliminates the patient's shift over the mattress and thus maintains a uniform position of the patient's body that is not bound to the position of the bed parts.

# 12 Scales Control (only version with scales)

Use iBoard Basic to control the scales.

## iBoard Basic




Fig. Scales Section (iBoard Basic) - display and keyboard

## 12.1 Preparation

- Install mattress and accessories to prepare bed before patient admission and using the scales.



**CAUTION!**  
**Incorrect use of scales due to incomplete preparation!**  
 ► Before each patient admission tare the scales.

## 12.2 Taring

Taring can be done in a range of 5kg to 249.5kg. Taring is used to set "0" on the display before placing the patient on the bed. Taring must be done with an unloaded bed with mattress, bed sheets, pillows and necessary accessories, without the patient. It is recommended to position Mattress support platform about 20 cm above the lowest horizontal position.

### To tare weight:

- Ensure that nothing and nobody touches the bed except you.
- Press and hold button  until value (6) starts to flash. Release button .
- Press button  again to confirm taring. "0" is shown on the display.

Place the patient on the bed.

### To cancel taring:

- Press button  while taring.

## 12.3 Displaying

Verification Scale Interval is 0.5 kg.

- Press button  to display value with actual scale interval 0,1 kg for 5s.

Unit of weight (**kg**) and decimal point are flashing on the display during this mode.

Field 6 shows actual weight value that remains displayed on the display.

If the bed is configured as EMR ready weight value automatically disappears after 1 minute and bed picture appears instead of it on the display. If it is needed, press button  to display weight value again.

## 12.4 Hold Mode

Hold Mode can be used only when scales are stabilized.

It allows adding or removing bed accessories and other items without changing the weight value.

### To activate Hold Mode:

- ▶ Wait until the scales are stabilized. The icon  or  will be illuminated when the scales are stabilized.
- ▶ Press button  until snowflake icon  appears on the display.
- ▶ Add or remove required accessories.

### To deactivate Hold Mode:

- ▶ After adding or removing accessories wait until the scales are stabilized (icon  or  is displayed).
- ▶ Press button .
- ▶ Display shows the original weight value.

### To deactivate Hold Mode without fixing the weight value:

- ▶ Press button .

## 12.5 Bed Overload

If load of the bed is over 254,5 kg:

- ▶ The „HIGH“ is shown on the display.

If load of the bed is over 260 kg:

- ▶ The „OVERLOAD“ is shown on the display.

**NOTE** If the bed is overloaded, it is impossible to position or manipulate the bed until overloading is removed.

**NOTE** Bed overloading always has higher priority than Hold Mode and Taring.

## 12.6 Bed Underload

If the bed is underloaded:

- ▶ Display shows the „LOW“.

## 12.7 Weighing in tilt

Accuracy is conditioned by the spirit level, which is located on the right head corner of the bed. If the bubble is in the highlighted circle then weighing is accurate.

## 12.8 Zeroing Scales

Zeroing is only possible in a range of  $\pm 5$ kg from factory zero.

Zeroing is used to reset weight on the display and set up user zero, which sets the maximum weight range of the weighing system. Zeroing must be done with an empty, unloaded bed, without the mattress and accessories. Zeroing is done after installation, weight verification or servicing.

### To zero scales:

- ▶ Position the bed about 20 cm above the lowest position and set the mattress support platform to the horizontal position. Ensure that nothing touches the bed except you.
- ▶ Press and hold button  until weight value starts to flash.
- ▶ Press button  to confirm zeroing.

“0” is shown on the display and an acoustic signal confirms zeroing.

### To cancel zeroing:

- ▶ Press button  while zeroing.



Fig. HOLD Mode (iBoard Basic)



Fig. High Load (iBoard Basic)



Fig. Bed Underload (iBoard Basic)

## 13 Bed Exit Monitoring (only version with scales)

Use iBoard Basic to control the Bed Exit Monitoring.

If the bed is equipped with SafetyMonitor (accessory), Bed Exit Monitoring function is included between the functions of the Safety-Monitor.

### iBoard Basic




Fig. Bed Exit Monitoring Section (iBoard Basic) - display and keyboard

### 13.1 Preparation

- Place a patient on the bed with suitable mattress.
- For the correct Bed Exit Monitoring in the Inner Zone patient's position in the middle of the bed is needed.

### 13.2 Activation of Bed Exit Monitoring

Bed Exit Monitoring is **OFF** and icon **OFF** is displayed by default.

To activate Bed Exit Monitoring:

- Press button .

Icon **ON** appears on the display.

When Bed Exit Monitoring is activated Inner Zone is set by default. Icon  therefore appears on the display.

**NOTE** Minimum patient weight for Bed Exit Monitoring is 35 kg.

### 13.3 Monitored Zone

To switch to the Outer Zone:

- Press button .

Icon  appears on the display.

To set the Inner Zone again:

- Press button .

Icon  appears on the display.

## 13.4 BED EXIT ALARM

Alarm is triggered when patient has left selected monitored zone or PAUSE period elapsed and patient is not in ordered position.

### To stop Alarm:

- Press button .

Bed Exit Monitoring is deactivated and icon **OFF** appears on the display.

The audible alarm is stopped.

### To pause Alarm:

- Press button  or button .

„Pause 15 min.“ appears on the display and 15 minute countdown timer runs.

The audible alarm is paused.



Fig. Visual signalisation of Bed Exit Alarm (iBoard Basic)

### 13.4.1 Alarm Volume

Maximum Alarm Volume Level is set by default.

It is possible to set Alarm Volume before and during triggered alarm.

#### To lower Alarm Volume Level:

- Press button .

Icon  with lower Alarm Volume Level appears on the display. Volume is lowered.

#### To return to Maximum Alarm Volume Level:

- Press button  after Minimum Alarm Volume Level has been reached.

Icon  with the 3 levels appears on the display.




Fig. Volume Icon

1. Minimum Volume
2. Moderate Volume
3. Maximum Volume

## 13.5 PAUSE

During PAUSE Mode Bed Exit Monitoring is temporarily interrupted and Alarms are not activated.

When Outer Zone monitoring is activated, PAUSE period is terminated when patient returns to the bed within 15 minutes of the PAUSE period.

### To PAUSE Bed Exit Monitoring:

- Press button  or button .

Icon **10** appears on the display with 15 minute countdown timer.

After PAUSE period elapsed and patient is in ordered position Bed Exit Monitoring is reactivated.

### To extend the PAUSE period:

- Press button  again to extend the countdown to 15 minutes period again.

### To terminate the PAUSE period:

- Press button .



Fig. PAUSE Countdown (with remaining minutes)

## 13.6 Deactivation of Bed Exit Monitoring

### To deactivate Bed Exit Monitoring:

- Press button .

Icon **OFF** appears on the display.





Fig. Bed Exit OFF (iBoard Basic)

## 14 Equipment

### 14.1 Accessory Rail with plastic hooks (optional)


Accessory Rail with 2 plastic hooks is intended for hanging accessories.  
It is located on the sides of bed.



*Fig. Accessory Rail with plastic hooks (on side)*

### 14.2 Linen Shelf (standard)

Linen Shelf is intended for storing the things corresponding to its dimensions (e.g. linens or sheets).  
The Attendant Control Panel can be stored in the Linen Shelf.  
The Linen Shelf is located at foot end under the Foot Board.  
Safe Working Load of the Linen Shelf is 15 kg.



*Fig. Linen Shelf (at foot end)*

## 14.3 Mobi-Lift® (optional)



### WARNING!

**Risk of injury due to slipping or falling when standing up!**

- ▶ Ensure that the Mobi-Lift handles are completely inserted in the sleeve fittings.
- ▶ Ensure that no bed linen is caught between the sleeve fitting and the support handle.

Mobi-Lift® is optional. It serves as a support handle to enhance the patient's safety when getting up. Mobi-Lift® is a support handle with a built-in Height Adjustment button. It allows the patient to raise and lower the Mattress support platform.

## Using the Mobi-Lift handles

### To adjust the Mobi-Lift handle:

- ▶ Lift the handle up towards the bed.
- ▶ Push the handle into the sleeve fitting as far as it will go.

### To adjust the height of the mattress support platform:

- ▶ Press button  on any control element.
- ▶ Press the button to adjust the bed height.




Fig. Mobi-Lift Handle

---

## 14.4 Brake Signal (optional)

If bed is equipped with Brake Signal and this bed is connected to the mains power, the Brake Signal sounds when the bed is not braked.

## 14.5 i-Brake® (optional)

It is possible to equip the bed with an automatic castor brake. The automatic castor brake prevents injuries of patients and staff due to an unbraked bed. The brakes are activated automatically 60 seconds after the bed is plugged in, and 60 seconds after they have been released if the bed is not being moved. It is possible to activate the brakes manually as well.

## 14.6 Retractable Fifth Castor (optional)

It is possible to equip the bed with Fifth Castor in the centre of undercarriage. The Fifth Castor helps to steer and manoeuvre the bed in long corridors and small rooms. Retracted Fifth Castor does not obstruct access to any devices under the undercarriage.

## 14.7 LINIS SafetyPort (optional)

LINIS SafetyPort is a medical device data system for capturing and transferring data from LINET beds into SafetyPort Dashboard and third party systems, including nurse calls, EHR and digital whiteboards. Data collection and evaluation takes place at one central location for all beds connected to the system simultaneously. The records are completely anonymous and the system does not work with any personally identifiable information. The customer can decide which data will be sent to the 3rd party system and adjust their sending period. LINIS SafetyPort is intended to be used to increase efficiency of healthcare personnel workflows by saving their time spent on documentation and eliminating errors. This is achieved by automated recording of different parameters of medical beds and their subsequent transfer to various hospital systems in HL7 format. Optional feature LINIS SafetyPort Dashboard is intended to save time the healthcare personnel spends on checking different beds at their workspace and to provide them with both near-real-time data and their aggregation to be able to check the history of provided care. LINIS SafetyPort may be used in various healthcare environments, including both intensive and non-intensive care units as well as units providing speciality care to a broad population of patients. The product is intended to be used by variety of healthcare personnel who have the cognitive skills to operate the product and are trained to use the product. LINIS SafetyPort is not an alarm system and the use of this product for this purpose means incorrect use.

## 14.8 i-Drive Power (optional)

### 14.8.1 i-Drive Power System - Basic Description

It is possible to equip the bed with the i-Drive Power wheel. The i-Drive Power helps hospital staff to drive the bed during patient transport with minimal manpower.

The i-Drive wheel is located in the center of the bed under the undercarriage. i-Drive Power is equipped with its own accumulator and charger and it is not dependent on the bed functions so, if discharged you can still use the bed functions. The bed is equipped with one i-Drive controller. i-Drive is oriented in straight direction of the bed.

### 14.8.2 Safety instruction for i-Drive Power

- ▶ Follow the instructions carefully.
- ▶ Ensure that the bed is operated exclusively by qualified staff.
- ▶ Make sure the siderails are raised up during the transport.
- ▶ Never use bed positioning buttons during transport.
- ▶ Never use Fast forward button when descending. The Fast forward button is recommended for use when ascending as it is more efficient.
- ▶ Special precaution need to be considered when reversing. Always keep distance from the bed and never use reverse button when descending or ascending.
- ▶ Do not use Free Drive to transport on a slope over 1 degree unless adequate personnel are available to manage safe bed transport.
- ▶ The driving down the slope that exceeds 6 degrees will require adequate contribution of a manpower.
- ▶ Never leave the bed with an activated i-Drive Power system without supervision of the trained staff.
- ▶ Always use the regular mechanical brake system to brake and stabilize the bed.
- ▶ Pay increased attention when driving the bed using i-Drive Power. Be aware of people and objects in close proximity and avoid collision with them by careful driving, especially by appropriate speed control.
- ▶ Make sure the bed is unplugged and bed brakes are released before using i-Drive Power.
- ▶ Push the emergency stop drive button if immediate movement interruption is needed (e.g. to avoid collision with other persons or objects).
- ▶ Retract the i-Drive Power wheel to the undercarriage when parking. This will prevent misuse when unbraking and braking the bed.
- ▶ The i-Drive Power electromagnetic brake is designed just for temporary bed stop and not for the permanent parking.
- ▶ Switch off the i-Drive Power accumulator prior to long-term storage or transport.
- ▶ Push the emergency retraction button under the chassis cover to retract the i-Drive Power wheel
- ▶ in case an of i- Drive Power system failure. This will enable moving the bed to a safe area manually without using i-Drive Power.
- ▶ Retract the i-Drive Power wheel to the undercarriage every time you intend to move the bed sideways.
- ▶ Pay attention to the LED accumulator status indicator and plan your drive using the i-Drive Power accordingly. Insufficient accumulator capacity can cause unexpected complications and risks during the drive.
- ▶ Always plug the bed in when you finish your drive in order to recharge the accumulator and keep your bed ready to go using the i-Drive Power.
- ▶ The i-Drive Power accumulator must be replaced every 2 years to maintain proper functions of the i-Drive Power.

### 14.8.3 Specifications of Use



#### **WARNING!**

##### **Risk of injury due to careless driving!**

- ▶ Always drive safely and carefully.
- ▶ Observe the path for any obstacles and avoid collisions.
- ▶ Ensure there are no people in your way.
- ▶ Manipulate with the bed carefully not to drive over any staff or patients.



#### **CAUTION!**

##### **Maximal clearance underneath the bed is 11,3 cm!**

- ▶ Observe the path for any obstacles and avoid collisions.

#### **Intended use:**

- ▶ bed transport (with or without patient) by the hospital staff

#### **Unintended use:**

- ▶ riding the bed
- ▶ other usage than described in instructions for use
- ▶ by other person than the trained staff

**NOTE** Each bed can transport only single patient at a time and cannot be used to transport other items (except bed accessories in secured position).

**NOTE** For information concerning uses other than those outlined in the “Specifications of Use” section above, please contact LINET®.

## 14.8.4 Manipulation



**CAUTION!**  
Damage to i-Drive Power main control panel cable due to wrong cable placement!

- Ensure that the main control panel connecting cable is placed correctly.



**CAUTION!**  
Material damage due to incorrect use!

- Do not hang anything on the main control panel and its cable!

1. Safety Sense (touch sensor)
2. Main Control Panel
3. Main Control Panel Cable – correct cable placement
4. Activation Panel




Fig. Position of Main Control Panel




Fig. Activation Panel (4)




Fig. Main control panel (2)

5. ON Button (i-Drive Wheel Activation Button)
6. OFF Button (i-Drive Wheel Retraction and Deactivation Button)
7. FAST FORWARD Button
8. STOP DRIVE Button
9. ON LED
10. Mains Power LED (bed connected to the mains or bed braked)
11. Fault LED
12. State of charge and Fault LED
13. FORWARD Button
14. REVERSE Button

**NOTE** The i-Drive Power controller cannot control the bed functions. Control the bed using the bed control elements.

**NOTE** The main control panel is enhanced with a touch sensor (1); your hand must always be in contact with the i-Drive Power control panel to use the functions. If released, the i-Drive Power will stop.

**NOTE** Raising and lowering of the i-Drive wheel is electrically controlled by the i-Drive activation panel.

## 14.8.5 i-Drive Power Activation/Deactivation

### To activate the i-Drive Power:

1. Check, if the mains switch of i-Drive Power is activated.
2. Press the Activation button  located on the Activation Panel. The i-Drive wheel will lower and the green indicator will flash.

### To deactivate the i-Drive Power:

1. Retract the i-Drive wheel using the button  located on the Activation Panel.
2. Deactivate the i-Drive using the mains switch.

### Emergency i-Drive Power wheel retraction:

1. Press any button  on the bed.
2. Deactivate the i-Drive Power using the mains switch.
3. Press the i-Drive Power Emergency Retraction Button situated on the bottom side of undercarriage under the label.

**NOTE** Use emergency retraction in case of accumulator discharge or drive malfunction to move the bed to a safe area manually without using i-Drive Power.





Fig. i-Drive Power Mains Switch with Label



1. i-Drive Power Mains Switch
2. i-Drive Power Emergency Retraction Button

### TO RETRACT THE I-DRIVE POWER WHEEL IN EMERGENCY:

1. PRESS ANY GO BUTTON ON THE BED
2. SWITCH OFF THE MAIN POWER SWITCH
3. HOLD THE EMERGENCY RETRACTION BUTTON UNDER THE CHASSIS COVER



Fig. i-Drive Power Emergency Retraction Button Label

## 14.8.6 Powered Drive



### CAUTION!

#### Damage to property due to incorrect transport and involuntary movement!

- ▶ Prior to transport, ensure that the bed is disconnected from the mains.
- ▶ Prior to transport, ensure that the auxiliary outlet plug (if available) is disconnected from the mains.
- ▶ Ensure that the castors are locked prior to assembly, disassembly and maintenance (e.g.: i-Drive Power maintenance).
- ▶ Ensure that the castors are locked when the bed is occupied.
- ▶ Hang the mains cable on the appropriate hook on the bed during transport.

1. Check, if the mains switch of i-Drive Power is activated.
2. Press the button  on the Activation Panel. The i-Drive wheel will lower and the LED  will flash.
3. Place your hand on the Safety Sense touch sensor (1) and push the button  or button  or button . Your hand must be placed on the Safety Sense sensor to use the i-Drive Power, if released, the i-Drive Power will stop.
4. The i-Drive motor is immediately stopped and the electric brake is activated after pressing the red button  when braking or in emergency.
5. i-Drive Power control system is automatically deactivated and the electric brake is activated if no i-Drive function is used for 3 minutes. This is signaled by the green LED  which is extinguished after 3 minutes.

**NOTE** i-Drive Power is not designed for ascending or descending a slope greater than 6° or longer than 20 m. The support of personnel is needed when ascending or descending with a full SWL.

**NOTE** When i-Drive wheel is lowered, it is not possible to move the bed sideways. Press the button  to retract the wheel, release the castors to the neutral position and then move the bed to any direction required.

## 14.8.7 Braking

1. Press and hold the button  to brake immediately.  
-or-
2. Press and hold the button  to brake slowly (Press the button  to brake when reversing).  
-or-
3. Release your hand from the touch sensor area and i-Drive Power will brake automatically.

**NOTE** Always brake the bed by using the castor control lever when the transport is finished or interrupted. The i-Drive electromagnetic brake is not designed to permanently brake the bed.

**NOTE** In a crisis situation (e.g. acceleration when driving down a steep slope) i-Drive dual braking prevents acceleration and slows down bed movement. However, it is not guaranteed the bed will stop by itself without personnel support (using button  and castor control lever).

**NOTE** When descending, it is possible to actively brake using the opposite direction button to slow.

## 14.8.8 Free Drive

The i-Drive motor is equipped with free drive, which is active after pressing the forwards (  or  ) or backwards (  ) buttons (until user holds the touch sensor area).

Free Drive is deactivated and the brake is activated when the direction of motion is changed. This is feature for lowering the risks when going to a slope.

## 14.8.9 Accumulator

**Accumulator charge status:**

1. While this indicator is flashing, the accumulator is critically discharged. (LED1)
2. 50% (LED2)
3. 75% (LED3)
4. 100% - the accumulator is charged (LED4)




Fig. Accumulator Charge Status

**To charge the accumulator:**

- Connect the bed main cable to mains power.
- i-Drive will be charged (with the accumulator discharged, the charging may take up to 9 hours).

**NOTE** Accumulator charge values are just informational. Accumulator life is reduced when the accumulator is allowed to discharge completely.

## 14.8.10 Fault Signalization

The system is protected against failure states, by stopping and braking the drive system, and respective signalization. The fault indicator flashing briefly and the accumulator indicator shows the fault state. Some defects are cleared automatically (e.g.: drive overheating). When drive or electronics is overheated, an short acoustic signal occurs before the drive is blocked.

| Error                                                 | LED1 | LED2 | LED3 | LED4 |
|-------------------------------------------------------|------|------|------|------|
| Drive overheated                                      | OFF  | OFF  | OFF  | ON   |
| Electronics overheated                                | OFF  | OFF  | ON   | OFF  |
| Brake error                                           | OFF  | OFF  | ON   | ON   |
| Retraction not completed                              | OFF  | ON   | OFF  | OFF  |
| Internal system error                                 | OFF  | ON   | OFF  | ON   |
| C losing of the Field-effect transistor is penetrated | OFF  | ON   | ON   | OFF  |
| Control circuit overheated                            | OFF  | ON   | ON   | ON   |
| Control circuit error                                 | ON   | OFF  | OFF  | OFF  |
| Activation button stuck                               | ON   | OFF  | OFF  | ON   |
| Retraction button stuck                               | ON   | OFF  | ON   | OFF  |
| Active button after start                             | ON   | OFF  | ON   | ON   |

## 14.8.11 Light Indicators

| Indicator                                                | Meaning                                                                                                                                                                                                                                                |
|----------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Go Indicator</b><br>► Constantly lit<br>► Flashing    | Hand is on touch sensor; drive wheel is ready for use.<br>Hand is not on touch sensor; i-Drive is not ready for use.                                                                                                                                   |
| <b>Fault Indicator</b><br>► Constantly lit<br>► Flashing | i-Drive cannot be activated (i-Drive wheel is not lowered, castor control lever is braked, bed is connected to the mains).<br>System is faulty (indicated on accumulator status indicator)<br>-or-<br>i-Drive control box heat protection is activated |

## 14.8.12 Technical Specifications

| Parameter                                     | Value            |
|-----------------------------------------------|------------------|
| i-Drive wheel diameter                        | 8,27 in.         |
| Max. fast forward speed (flat ground, loaded) | 4,43 Km/h (±15%) |
| Max. forward speed (flat ground, loaded)      | 2,16 Km/h (±15%) |
| Max. reverse speed (flat ground, loaded)      | 2,16 Km/h (±15%) |
| Max. angle of ascent                          | 6°               |
| Noise level (when retracting the drive wheel) | 65 dB            |

## 14.8.13 Electrical specification

| Parameter                | Value                                                                                                      |
|--------------------------|------------------------------------------------------------------------------------------------------------|
| Input Voltage, Frequency | 230 V AC, 50/60 Hz<br>127 V AC, 50/60 Hz<br>120 V AC, 50/60 Hz<br>110 V AC, 50/60 Hz<br>100 V AC, 50/60 Hz |
| Accumulator Voltage      | 36 V DC, Capacity: 12 Ah                                                                                   |
| Maximum Power Input      | 300 W                                                                                                      |
| Fuse                     |                                                                                                            |
| Version 230 V            | 2 x T1,6A L 250V                                                                                           |
| Version 127 V            | 2 x T3,15A L 250V                                                                                          |
| Version 120 V            | 2 x T3,15A L 250V                                                                                          |
| Version 110 V            | 2 x T3,15A L 250V                                                                                          |
| Version 100 V            | 2 x T3,15A L 250V                                                                                          |

## 14.8.14 i-Drive Power Maintenance

Periodical maintenance of the i-Drive Power must be done by qualified service technician or authorized service organization at least once a year. To continue maintenance please see chapter Maintenance.

### Service technician must check the following:

- ▶ accumulator status and eventual replacement of the accumulator (after maximum of three years of duty)
- ▶ gas spring – replace if necessary (after maximum of three years of duty)
- ▶ i-Drive Power wheel – replace if necessary
- ▶ lifting mechanism – grease if necessary
- ▶ cables, control elements – replace if necessary
- ▶ i-Drive Power function

## 14.9 Safestop (optional)

Safestop prevents user of the bed from injuries due to crushing by the lowered Mattress support platform.

When obstacle occurs on the undercarriage and Mattress support platform is going down, the motion is automatically stopped.

iBoard Standard display shows **SAFE STOP** +  and beeping is performed.

## 14.10 X-Ray Examination (optional)



### WARNING!

#### Respect maximum dimensions of x-ray cassettes!

- ▶ Maximum dimensions of any x-ray cassette for X-Ray Cassette Holder are 46,5 cm x 39 cm x 1,8 cm!



### WARNING!

#### Prevent x-ray cassettes from being damaged!

- ▶ Do not leave any x-ray cassette in the X-Ray Cassette Holder if X-Ray Examination should not be performed!



### CAUTION!

#### Prevent x-ray images from being devalued!

- ▶ Take x-ray images when the bed is braked and no part of the bed is forced to move!
- ▶ Follow the instructions in this instructions for use how to use the X-Ray Cassette Holder!

The Backrest of the bed consists optionally of HPL and is x-ray translucent. The bed is equipped with an x-ray cassette holder inserted under the Backrest left side. This design allows taking x-ray images of the patient's lungs without moving the patient manually.

## Necessary Steps before the Examination

- ▶ Make sure that patient is in centre of bed.
- ▶ Make sure that backrest is in lowest position and siderails are raised up.
- ▶ Pull out x-ray cassette holder.
- ▶ Insert x-ray cassette in the horizontal position.
- ▶ Insert back x-ray cassette holder with x-ray cassette so that the cassette centre indicator is exactly under the edge of the mattress support platform.
- ▶ Correct position of x-ray cassette holder using the tooth mechanism so that the upper edge of the x-ray cassette is exactly under the patient's shoulder line. For the correct orientation use the scale on the label. Indicate the position of the patient's shoulder line using the numbers on the scale. Move the x-ray cassette holder in such position so that the centre of the handrail is on the respective scale number.
- ▶ Adjust parameters of the x-ray device and do the image.

**NOTE** This procedure is above all suitable for patients who cannot be moved due to critical conditions (e. g. internal bleeding) or unstable patients.





Fig. X-Ray Cassette Holder


## 14.11 Undercarriage Cover




**CAUTION!**

**Risk of material damage due to objects on the undercarriage cover!**

- Do not place objects on the undercarriage cover!



*Fig. 2-part undercarriage cover (standard)*



*Fig. 1-part undercarriage cover (optional)*




Fig. 3-part undercarriage cover (standard)

## 14.12 USB Connector (optional)



### WARNING!

#### Risk of injury due to incorrect use!

- Ensure accessory plugged in USB connector is in pristine condition!
- User of the bed is responsible for the fact that this requirement is met.



### CAUTION!

#### Risk of material damage due to incorrect use!

- Do not plug heating element into USB connector!
- User of the bed is responsible for the fact that this requirement is met.



### CAUTION!

#### Risk of material damage due to incorrect cleaning!

- No liquid should get to the USB port!

USB Connector situated on the both sides of Backrest is intended for charging mobile phones and tablets.

**NOTE** Maximum current for this device is 2 A.




Fig. Fig. USB label with electric specifications




Fig. USB Connector on the side of Backrest

# 15 Mattress

Eleganza 4 bed is designed for mattresses from LINET portfolio.



**CAUTION!**

**Incompatibility with bed due to incorrect mattress dimensions!**

- ▶ Check maximum approved mattress dimensions (chapter Technical Specification of Eleganza 4).

The manufacturer recommends the use of the following mattresses on the Eleganza 4 bed:

**PASSIVE MATTRESSES**


- CliniCare 10
- CliniCare 20
- CliniCare 30

**ACTIVE MATTRESSES**

- CliniCare 100 HF
- Air2Care
- Virtuoso

## 15.1 Installation of Passive Mattress

The passive mattresses intended for Eleganza 4 bed are equipped with straps with buckles (1) to fix mattress on the Mattress Support Platform.



*Fig. Bottom of Passive Mattress*

### 15.1.1 Straps with side release buckles

**To fix mattress on the Mattress Support Platform:**

- ▶ Run three straps through the three corresponding holes in the covers of Mattress Support Platform.
- ▶ Run these three straps under the bars of the Mattress Support Platform.
- ▶ Lock the three side release buckles by connecting their male and female parts together.

**To remove mattress from the Mattress Support Platform:**

- ▶ Release the three buckles by pressing them from both sides and by disconnecting their male and female parts.
- ▶ Pull these three straps out of the Mattress Support Platform.
- ▶ Remove mattress from the Mattress Support Platform.




Fig. Fixation of the passive mattress with straps on the mattress support platform of Eleganza 4 bed

## 15.2 Installation of Active Mattress



### WARNING!

Follow instructions for use of a compatible active mattress carefully!



### CAUTION!

Risk of material damage due to an incorrect fixation of compatible active mattress on the mattress support platform!

- ▶ Adjust the bed to maximum Cardiac Chair Position before fixing all the straps of the inflated mattress to the mattress support platform!

### Installation instructions:

- ▶ Remove any existing mattress.
- ▶ Observe mattress dimensions and its orientation before putting it on the Mattress support platform.
- ▶ Place SCU on the foot board of the bed or on the floor.


### 15.2.1 Straps with side release buckles

#### To fix mattress on the Mattress Support Platform:

- ▶ Inflate the mattress intended to be fixed on the mattress support platform.
- ▶ Adjust the bed to the maximum Cardiac Chair Position.
- ▶ Run four straps through the four corresponding holes in the covers of Mattress Support Platform and two straps through the accessory holders under the calfrest cover.
- ▶ Lock all the side release buckles by connecting their male and female parts together.

#### To remove mattress from the Mattress Support Platform:

- ▶ Release all the buckles by pressing them from both sides and by disconnecting their male and female parts.
- ▶ Pull all the straps out of the Mattress Support Platform.
- ▶ Remove mattress from the Mattress Support Platform.



*Fig. Fixation of active mattress with straps on the mattress support platform of Eleganza 4 bed*

## 16 Accessories



### **WARNING!**

#### **Risk of injury due to incompatible accessories!**

- Use exclusively original accessories from the manufacturer.
- The manufacturer is not responsible for the use of unapproved accessories.



### **WARNING!**

#### **Risk of injury due to damaged accessories!**

- Use exclusively accessories in perfect condition.



### **WARNING!**

**Avoid collisions between accessories and bed parts before use of the accessories!**




### **WARNING!**

#### **Risk of injury or material damage due to incorrect use!**

- Compatible accessories manufactured by different manufacturers have their own instructions for use. It is necessary to read instructions for use of a compatible accessory with instructions for use of the compatible LINET product to respect especially technical parameters, warning notifications, cleaning and maintenance instructions of LINET products and their compatible accessories!

| Compatible Accessories                                                | Identification Numbers                                               |
|-----------------------------------------------------------------------|----------------------------------------------------------------------|
| <b>Lifting pole</b>                                                   | 11011450A0000, 11011450B0000, 1101145000000                          |
| with Triangular holder                                                | 4ROTGERSG700-2, 4ROTGNR100GR-2                                       |
| with Infusion holder                                                  | 16010700B0000, 16010700A0000, 4MAPL00N1001, 1101030000000            |
| <b>Telescopic Infusion Stand</b>                                      | 4PV17744500L, 4PV340290000, 4PV322717000, 4PV17604800L, 4PV17587200L |
| with Fixation holder                                                  | 11028700B0000                                                        |
| with Infusion bottle basket                                           | 4DR426101                                                            |
| <b>Collapsible Infusion Stand</b>                                     | 4MA45MNA0006                                                         |
| <b>Mobile Infusion Stand</b>                                          | 4ZZ426100                                                            |
| <b>Writing shelf</b>                                                  | 11023230B0000                                                        |
| <b>Monitor shelf</b>                                                  | 11026700B0000                                                        |
| <b>Utility shelf</b>                                                  | 11025800B0000                                                        |
| <b>Crutches holder</b>                                                | 1102070000000                                                        |
| <b>Name holder</b>                                                    | 4GT006000000                                                         |
| <b>Horizontal oxygen bottle holder</b>                                | 4MAS6015553                                                          |
| <b>Vertical oxygen bottle holder</b>                                  | 4MAR2010PC007-1                                                      |
| with Adaptor                                                          | 11029700A0000                                                        |
| <b>Siderail accessory holder</b>                                      | 11028700B1000                                                        |
| <b>Ventilation circuit holder</b>                                     | 4MAR2016PC001                                                        |
| <b>Urinary bag holder</b>                                             | 11023820B0000                                                        |
| <b>Positioning cushion</b>                                            | 4SET17500001                                                         |
| <b>Traction frame E</b>                                               | 4MAEXE500000                                                         |
| <b>m-Panel</b>                                                        | 1MPAS6018146                                                         |
| <b>Pair of Protectors</b>                                             | 11005900B0ELE4                                                       |
| <b>Pushing Handles</b>                                                | 11022200A0000                                                        |
| <b>Pushing Handles with Holder of the i-Drive Power Control Panel</b> | 11022200A0X01                                                        |

1. Lifting pole
2. Infusion stand
3. Oxygen bottle holder (vertical version)
4. Ventilation circuit holder



*Fig. Accessories*

## 16.1 Lifting Pole



### WARNING!

Risk of injury or material damage due to use of the Lifting Pole 11011410B0000 or 11011410A0000 or 1101141000000 with Eleganza 4 bed!

- Avoid collisions between the Lifting Pole 11011410B0000 or 11011410A0000 or 1101141000000 and the head board of the bed in the Trendelenburg position!

### To ensure safe use of the lifting pole:

- Never exceed the maximum load of 75 kg.
- Never use the lifting pole for rehabilitation exercises.
- To prevent the bed from tipping over, ensure that the lifting pole does not project out from the bed.
- Replace plastic handle every 4 years.

### To install the lifting pole:

- Insert lifting pole in corresponding sleeve fitting on accessory adapter at head end.
- Ensure that safety pin locks into place.
- Attach a plastic grab handle with an adjustable strap to the lifting pole.

**NOTE** The date of manufacture is marked on the grab handle.




Fig. Places for lifting pole (sleeve fittings on accessory adapter)




Fig. Lifting pole

## 16.2 Infusion Stand



### **WARNING!**

#### **Risk of injury due to use of incorrect accessories or because of incorrect use!**

Infusion Stands must only be used for their intended use. Always read the instructions for use!

- ▶ Only mount an infusion pump to the lower (wider) telescopic section of an infusion stand above the head/foot board.
- ▶ Never mount an infusion pump to the upper (thinner) telescopic section of an infusion stand.
- ▶ Ensure the infusion pump will not collide with any movable parts of the bed (especially Backrest part) or with the patient. This must be verified during installation.
- ▶ Do not over tighten the infusion pump clamps during fitment. Over tightening may damage the infusion stand.
- ▶ Infusion pump can be only used if the infusion stand is fitted in the accessory holder socket in the head end on the under carriage of the bed.
- ▶ Do not use the infusion stand as driving/pushing device during the bed transport.



### **CAUTION!**


#### **Risk of collision with oxygen bottle holder on the bed end due to incompatibility!**

- ▶ Use the foldable infusion stand with adapter to avoid the collision.


Infusion stands are intended to provide a suitable support for the attachment of the infusion pumps/syringe pumps and suspension of the infusion bags or bottles.

Infusion stands can be fitted to the head and foot end of the bed by either fitting into the IV/Infusion sockets mounted on the bed or using alternative accessory holder socket in the head end on the undercarriage of the bed.

- ▶ Use exclusively infusion stands with hooks for hanging IV bags or baskets for intravenous solutions.
- ▶ Ensure the infusion stand individual hook 2kg maximum Safe Working Load is not exceeded.
- ▶ Ensure the infusion stand maximum Safe Working Load is not exceeded.
- ▶ Follow the actual price list for information about types of infusion stands.



*Fig. Places for infusion stand (sleeve fittings on accessory adapter)*



*Fig. Infusion stand*

## 16.3 Oxygen Bottle Holder



### WARNING!

**Risk of injury with oxygen bottle holder due to incorrect use or due to careless driving!**

- ▶ Ensure the oxygen bottle holder is correctly fitted in correct position.
- ▶ It is necessary to place oxygen bottle holder (with or without oxygen bottle) before transport to secure transport position.
- ▶ Be aware of people or objects in close proximity when driving or manipulating the bed equipped with oxygen bottle holder.
- ▶ Secure the oxygen bottles against falling or involuntary movement with rubber strap.
- ▶ Place the oxygen bottle holder on the bed by instructions in the following text.
- ▶ Ensure the oxygen bottle valve will not get damaged by careless or incorrect manipulation or placement.

The oxygen bottle holders are suitable for transporting oxygen bottles with a weight of up to 15 kg and a volume of 5 litres.

### Horizontal version (Head End)

- ▶ Put oxygen bottle holder on crossbar on the head end of the bed.




Fig. Oxygen Bottle Holder (on the head end)

### Vertical version (with adapter)

On following pictures there are 4 positions of the oxygen bottle holder with adapter.

- ▶ Put holder on sleeve fittings in multifunctional accessory adapter on head end.
- ▶ Ensure the locking pin of vertical oxygen bottle holder is locked in sleeve fitting.




Fig. Oxygen Bottle Holder in the adapter (position 1)




Fig. Oxygen Bottle Holder in the adapter (position 2)




Fig. Oxygen Bottle Holder in the adapter (position 3)




Fig. Oxygen Bottle Holder in the adapter (position 4)

## 16.4 Ventilation Circuit Holder

The ventilation circuit holder prevents extubation of the patient connected to the ventilator.

- ▶ Always use LINET ® ventilation circuit holder to prevent extubation during any procedures.

### Applying ventilation circuit holder:

- ▶ Put ventilation circuit holder in hole on right or left side of the Backrest frame.
- ▶ Fasten ventilation circuit holder with wing screw provided.
- ▶ Put intubation tube through plastic head of ventilation circuit holder.
- ▶ Tilt mattress support platform left and right by 15° to check if intubation tube is fastened securely. The fastening is secure if no parts of the ventilation circuit are disconnected.





Fig. Ventilation Circuit Holder

## 16.5 Writing Shelf

The Writing Shelf is intended for writing of nursing staff.  
It is placed in the handles of the foot board (as on the picture).




*Fig. Writing Shelf*

## 16.6 Monitor Shelf

The Monitor Shelf is suitable for transporting monitors with a weight of up to 15 kg.

### Installing the Monitor Shelf:

- ▶ Place Monitor Shelf on the foot board (as on the picture).
- ▶ Fix the monitor with safety belts in order to avoid any damage during transport.



*Fig. Monitor Shelf*

## 16.7 Urinary Bag Holder

Urinary Bag Holders are available on both sides of the bed at Backrest area.




Fig. Urinary Bag Holder

## 16.8 Protectors



### **WARNING!**

#### **Risk of injury due to the patient falling off the bed!**

- ▶ Ensure that the Protector is installed securely.
- ▶ Always check that the siderails are properly locked in the „up“ position.
- ▶ Make sure the fall risk assessment was done properly before use of the Protectors.
- ▶ For use of the Eleganza 4 bed with Protectors, it is recommended to disable the Backrest and Thighrest Adjustment functions of the Patient Control Panels located on the inner sides of the foot siderails via Attendant Control Panel or Nurse Control Panel to reduce the risk of patient fall to the minimum level.



### **CAUTION!**

**Right Protector and left Protector are not interchangeable! Ensure the both Protectors are properly installed on the bed!**



### **CAUTION!**

**Follow the wash symbols on the label of Protector cover to clean the Protector cover correctly!**

## Intended Use

Two Protectors are an optional accessory for the Eleganza 4 bed. The Protectors can only be installed on the Eleganza 4 bed ready for Protectors. This bed configuration includes the casings for Protectors on both sides of the bed at foot end. Without this bed configuration, an intervention of a qualified and trained service technician is necessary to enable the use of the Eleganza 4 bed with Protectors. The main purpose of the Protectors is to reduce the risk of fall at very risky patients (confused restless patients). The use of the Protectors for the other patients is not necessary. The Protectors can be used with standard or extended mattress support platform.




Fig. Protector placement at the right foot end of the Eleganza 4 bed


## Placement

### Attach the Protector to the bed as follows:


- ▶ Remove the Foot Board (1).
- ▶ Extend the Mattress Support Platform as needed (2).
- ▶ Lift the Calfrest up (3).
- ▶ Turn the safety catch (4) to ensure the casing for the Protector is opened.
- ▶ Insert the Protector to the casing at foot end of the bed with regard to the required Mattress Support Platform extension (5).
- ▶ Turn the safety catch back to ensure the casing for the Protector is closed in required position.
- ▶ Lower the Calfrest down.
- ▶ Reduce the Mattress Support Platform extension as needed.
- ▶ Insert the Foot Board back to the foot end of the bed.

### Remove the Protector from the bed as follows:

- ▶ Remove the Foot Board (1).
- ▶ Extend the Mattress Support Platform as needed (2).
- ▶ Lift the Calfrest up (3).
- ▶ Turn the safety catch (4) to ensure the casing for the Protector is opened.
- ▶ Remove the Protector from the casing at foot end of the bed.
- ▶ Lower the Calfrest down.
- ▶ Reduce the Mattress Support Platform extension as needed.
- ▶ Insert the Foot Board back to the foot end of the bed.



*Fig. Position of the right safety catch for fixing the right Protector*



*Fig. Protector correctly installed at the right foot end of the Eleganza 4 bed*

## 16.9 SafetyMonitor



### WARNING!

Read the instructions for use for SafetyMonitor carefully before use of the Eleganza 4 bed with SafetyMonitor!

### Intended Use

SafetyMonitor is a system that monitors statuses of Eleganza 4 bed: brake status (braked castors/unbraked castors), siderail status (siderails up/siderail down), bed height (bed in the lowest position/bed not in the lowest position), backrest angle (backrest in more than 30°/backrest in less than 30°), patient presence (patient on the bed/patient not on the bed) and location of the bed (where the bed is situated). Secure statuses are: braked bed, siderails up and locked, bed in the lowest position, backrest in more than 30° and patient on the bed. The system triggers alerts (notice signals for insecure statuses) and alarms (alarm signals for the absence of a patient on the bed). Alerts and alarms are automatically passed to the hospital information system and displayed on a screen in the nurse station and on a smartphone/tablet. Information is transmitted via LAN or Wi-Fi connection. Alarm (Bed Exit Alarm) sounds just from the bed. In this way, medical staff can be informed in a timely manner of any safety risk while saving administrative time.

### System components

SafetyMonitor system consists of installed server, secure intranet infrastructure (Wi-Fi or LAN), screen in the nurse station (PC or tablet or smartphone), Tag on the wall, parked Eleganza 4 EMR ready bed with iBoard Basic, Integration Module, Localisation Receiver, LAN connector and LAN cable.

### iBoard Basic (SafetyMonitor Section)




Fig. SafetyMonitor Section (iBoard Basic) - Display and Keyboard

1. ON button (RUN)
2. Inner Zone button (Bed Exit Monitoring)
3. Outer Zone button (Bed Exit Monitoring)
4. OFF button
5. Volume button (3 levels of sound pressure)
6. PAUSE button
7. Bed Exit Monitoring Activated (Inner Zone Monitoring)
8. Patient icon (patient is on the bed)
9. Volume indicator (3 levels of sound pressure)
10. Server connection icon (arrows only – connected, arrows with cross – disconnected)
11. Bed icon with bed statuses
12. SafetyMonitor system status indicator (ON/OFF/PAUSE)




Fig. Bed icon with bed statuses  
(orange - alerts/insecure statuses)

1. siderail status (orange - siderail down)
2. backrest angle (orange - backrest in less than 30°)
3. bed height (orange - bed not in the lowest position)
4. brake status (orange - unbraked castors)

## 16.10 Traction Frame E



### WARNING!

#### Risk of injury due to incorrect use!

- ▶ There must be no patient on the bed when installing the Traction Frame E!
- ▶ Avoid collisions between the Traction Frame E and the bed (Foot Board, Backrest) during bed positioning!
- ▶ Avoid collisions between the Traction Frame E and accessories!
- ▶ Cross the thresholds with caution during transporting the bed with installed Traction Frame E!
- ▶ It is allowed to transport a patient on the bed with Traction Frame E just in emergency cases and with caution!
- ▶ Respect the Safe Working Load of the bed, of the Traction Frame E and of its hooks and pulleys!
- ▶ Remove the Traction Frame E from the bed if it is not needed for treatment!

### Intended Use

Traction Frame E is a supporting construction intended for fixation, traction and relief of limbs, spine and pelvis.  
Traction Frame E is intended for orthopaedic department, surgery department, traumatologic department and for ICU.

### Placement

Traction Frame E is inserted to the holes in the Accessory Adaptor at head end and to the holes in the Traction Frame Holder at foot end.




Fig. Eleganza 4 with Traction Frame E (side view)




Fig. Safe Working Load of the hooks  
(Infusion Holder)




Fig. Safe Working Load of the pulley




Fig. Eleganza 4 with Traction Frame E (top view)




Fig. Positions to insert Traction Frame E (undercarriage of Eleganza 4)

## 16.11 Pushing Handles



### **WARNING!**

**Pushing handles are not compatible with i-Drive Power system!**

- ▶ Do not use the pushing handles on the bed that is equipped with i-Drive Power system!




### **CAUTION!**

**Follow these instructions during placement of the pushing handles on the head end crossbar:**


- ▶ Placement of the pushing handles must be performed by hospital technician according to these instructions for use!
- ▶ Pushing handles can only be placed on the head end crossbar!
- ▶ Positions of the pushing handles are shown on the picture below!
- ▶ Distance between bushing holder and the bar that is at right angle to the head end crossbar is 12 mm.
- ▶ Nuts of bolts must be oriented inwards!

**Pair of Pushing Handles is intended for bed transport.**  
**Pushing Handles are not compatible with head board.**  
**Pushing Handles are located on the head end crossbar.**



**12 mm**

*Fig. Distance between bushing holder and the bar that is at right angle to the head end crossbar*



*Fig. Placement of pushing handles on the head end crossbar (positions of nuts)*

## Removal of the Pushing Handles from the bushings:

- Pull both Pushing Handles out of the fixed bushings on the head end crossbar.

## Insertion of the Pushing Handles to their fixed bushings:

- Insert both Pushing Handles to their fixed bushings on the head end crossbar.




Fig. Pushing Handles on the Eleganza 4 bed

## 16.12 Pushing Handles with Holder of the i-Drive Power Control Panel



### CAUTION!

Follow these instructions during placement of the Pushing Handles with Holder of the i-Drive Power Control Panel on the head end crossbar:


- Placement of the pushing handles must be performed by hospital technician according to these instructions for use!
- Pushing handles can only be placed on the bed head end crossbar!
- Positions of the pushing handles are shown on the picture below!
- Distance between bushing holder and the bar that is at right angle to the head end crossbar is 12 mm.
- Nuts of bolts must be oriented inwards!

Pair of Pushing Handles with Holder of the i-Drive Power Control Panel is intended for bed transport.

Pushing Handles with Holder of the i-Drive Power Control Panel are not compatible with bed head board.


Pushing Handles with Holder of the i-Drive Power Control Panel are located on the head end crossbar.

Pushing Handles with Holder of the i-Drive Power Control Panel are compatible with i-Drive Power system. i-Drive Power Control Panel is intended to be placed on the Holder of the i-Drive Power Control Panel.



12 mm

Fig. Distance between bushing holder and the bar that is at right angle to the head end crossbar




*Fig. Placement of pushing handles on the head end crossbar (positions of nuts)*

**Removal of the Pushing Handles from the bushings:**

- Pull both Pushing Handles out of the fixed bushings on the head end crossbar.

**Insertion of the Pushing Handles to their fixed bushings:**

- Insert both Pushing Handles to their fixed bushings on the head end crossbar.



*Fig. Pushing Handles with Holder of the i-Drive Power Control Panel on the Eleganza 4 bed*

## 16.13 IDock (size S and size M)



### WARNING!

#### Risk of injury or material damage due to incorrect use!

- ▶ Read Instructions for use of the IDock docking-cart before use of this product!
- ▶ Ensure the safe working load of the IDock docking-cart is not exceeded!



### CAUTION!

#### Risk of collision due to incompatibility!

- ▶ Ensure the selected type of IDock is compatible with the bed undercarriage!
- ▶ IDock is not compatible with any active mattress control unit placed on the Foot Board!



### CAUTION!

#### Risk of injury or material damage due to incorrect installation!

- ▶ The hospital personnel has to ensure the secure mounting of devices on the IDock and of the IDock on the bed!

IDock is intended to transport medical devices and oxygen bottles.

IDock is compatible with Multicare X bed, Multicare bed, Eleganza 5 bed and Eleganza 4 bed.

IDock size S and IDock size M are 2 types of the IDock docking-cart.

IDock must be connected with Foot Board via 2 hooks. If IDock is equipped with Transport Bracket (WDE2903860000) installed on the Foot Board of Multicare X bed or Multicare bed, IDock should be connected with the Transport Bracket via 2 hooks instead.

## Installation of the Transport Bracket (only for Multicare X and Multicare)



### WARNING!

Ensure the Transport Bracket is correctly installed on the bed foot end crossbar before connecting IDock to the Transport Bracket!





Fig. Positions for IDock hooks on the Transport Bracket



Fig. IDock docking-cart



*Fig. Installation of the Transport Bracket on the bed foot end crossbar*

## 17 Cleaning/Disinfection



### WARNING!

**Risk of injury due to accidental bed movement!**

- ▶ Always disable the function buttons when cleaning between the undercarriage and mattress support platform.



### CAUTION!

**Material damage due to incorrect cleaning/disinfection!**

- ▶ Do not use washing machines.
- ▶ Do not use pressure or steam cleaners.
- ▶ Follow the instructions and observe the dosages recommended by the manufacturer.
- ▶ Ensure that disinfectants are selected and applied exclusively by qualified hygiene experts.
- ▶ **Check if used cleaning agents and disinfectants are compatible with materials that the product consists of! For information see the following table.**

| BED COMPONENTS THAT ARE INTENDED TO BE CLEANED                                             | MATERIALS (SURFACES OF THE MENTIONED BED COMPONENTS)                                                                       |                                                                  |
|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------|
| Do not clean what is not mentioned in this column!                                         | Competent user is responsible for check if used cleaning agents and disinfectants are compatible with mentioned materials! |                                                                  |
| Head board and foot board                                                                  | Polypropylene (PP)                                                                                                         |                                                                  |
| Head siderails and foot siderails                                                          | Polypropylene (PP)                                                                                                         |                                                                  |
| Mattress support platform covers (Backrest)                                                | Polypropylene (PP)                                                                                                         | version with x-ray cassette holder: High Pressure Laminate (HPL) |
| Mattress support platform covers (Thighrest, Calfrest)                                     | Polypropylene (PP)                                                                                                         |                                                                  |
| Seat section                                                                               | Lacquered steel                                                                                                            |                                                                  |
| Castors                                                                                    | Polyurethane (PUR) + Polypropylene (PP)                                                                                    |                                                                  |
| Castor control levers                                                                      | Polyamide (PA6) + Lacquered steel                                                                                          |                                                                  |
| Frame of the mattress support platform                                                     | Lacquered steel                                                                                                            |                                                                  |
| Columns                                                                                    | Oxidized aluminium alloy                                                                                                   |                                                                  |
| Undercarriage cover                                                                        | Acrylonitrile butadiene styrene (ABS)                                                                                      |                                                                  |
| Corner covers                                                                              | Polypropylene (PP)                                                                                                         |                                                                  |
| Corner bumpers                                                                             | Polypropylene (PP)                                                                                                         |                                                                  |
| Keyboards (Attendant Control Panel, Handset, control elements integrated in the siderails) | Polyethylene terephthalate (PET)                                                                                           |                                                                  |
| CPR levers                                                                                 | Lacquered steel                                                                                                            |                                                                  |
| Labels                                                                                     | Polyethylene terephthalate (PET)                                                                                           |                                                                  |
| Accessory rail                                                                             | Polyoxymethylene (POM) + Lacquered steel                                                                                   |                                                                  |
| Actuators                                                                                  | Polyamide (PA6) + Aluminium (Al)                                                                                           |                                                                  |
| Decors (head board, foot board, head siderails, foot siderails)                            | Acrylonitrile butadiene styrene (ABS)                                                                                      |                                                                  |
| Mobi-Lift® handles                                                                         | Polyamide (PA66) + Lacquered steel                                                                                         |                                                                  |
| Foot controllers                                                                           | Acrylonitrile butadiene styrene (ABS) + rubber + Lacquered steel                                                           |                                                                  |

**For safe and gentle cleaning:**

- ▶ Do not use any strong acids or bases (optimum pH range 6 - 8).
- ▶ Exclusively use detergents that are suitable for cleaning medical equipment.
- ▶ Do not use abrasive powders, steel wool, or other materials and cleaning agents that might damage the mattress replacement system.
- ▶ Never use any corrosive or caustic detergents.
- ▶ Never use detergents that deposit calcium carbonate.
- ▶ Never use detergents with solvents that might affect the structure and consistency of the plastics (benzene, toluene, acetone, etc.).
- ▶ Clean electrical components carefully and allow them to dry completely.
- ▶ Do not immerse control unit in water or steam-clean it.
- ▶ Observe local directives regarding infection control.
- ▶ Make sure any cleaning agent used is approved by:
  - the facility in which the mattress replacement system is to be used.
  - by competent authority of the country in which the mattress replacement system is to be used.

## 17.1 Cleaning (Eleganza 4)

**Prepare for cleaning as follows:**

- ▶ Put the mattress support platform in the highest position.
- ▶ Adjust the back and thigh rests so that the reverse sides are accessible.
- ▶ Disable the function buttons on the control elements using the Attendant Control Panel.
- ▶ Disable the foot controls using the Attendant Control Panel.
- ▶ Disconnect the bed from the mains.
- ▶ Move the bed to the location where it will be cleaned.
- ▶ Lock the brakes on the bed.

### 17.1.1 Daily Cleaning

**Clean the following bed parts:**

- All control elements for adjusting the bed
- All handles
  - CPR release handle
- Bed ends
- Siderails (in highest position)
- Freely accessible mattress surface
- Mobi-Lift®
- Accessory rails

### 17.1.2 Cleaning before Changing Patients

**Clean the following bed parts:**

- All control elements for adjusting the bed
- All handles
- CPR release handle
- Bed ends
- Siderails (in highest position)
- Freely accessible mattress surface
- Mobi-Lift®
- Accessory rails
- All plastic mattress support platform covers
- Plastic undercarriage covers
- Telescopic columns
- Mattress on all sides
- Freely accessible metal parts of mattress support platform
- Cable ducts
- Lifting pole sleeve fitting
- Infusion stand sleeve fitting
- Bumpers
- Castors
- Brakes

## 17.1.3 Complete Cleaning and Disinfection

Clean the following bed parts:

- All control elements for adjusting the bed
- All handles
- CPR release handle
- Bed ends
- Siderails (in highest position)
- Freely accessible mattress surface
- Mobi-Lift®
- Accessory rails
- All plastic mattress support platform covers
- Plastic undercarriage covers
- Telescopic columns
- Mattress on all sides
- Freely accessible metal parts of mattress support platform
- Cable ducts
- Lifting pole sleeve fitting
- Infusion stand sleeve fitting
- Bumpers
- Castors
- Brakes
- Interior parts (accessible after removing mattress support platform covers)

## 18 Troubleshooting



**DANGER!**

**Risk of mortal injury due to electric shock!**

- If a fault occurs, have the electric motor, power box or other electrical parts repaired by qualified personnel exclusively.
- Do not open the protective covers of the electric motor or the power box.

| Error/Fault                                              | Cause                                                                   | Solution                                                      |
|----------------------------------------------------------|-------------------------------------------------------------------------|---------------------------------------------------------------|
| Adjusting with position buttons not possible             | GO Button was not pressed                                               | Press the GO button.                                          |
|                                                          | Function disabled on Attendant Control Panel                            | Enable disabled function.                                     |
|                                                          | Actuators have no power<br>Defective actuators<br>Defective accumulator | Check the mains connection.<br>Notify the service department. |
|                                                          | Mains Plug inserted incorrectly                                         | Insert the Mains Plug correctly.                              |
|                                                          | Faulty Power Source                                                     | Notify the service department.                                |
|                                                          | Faulty Control Element                                                  | Notify the service department.                                |
| Faulty Mattress support platform Height/Tilt Adjustment  | Obstacle on the undercarriage cover                                     | Remove the obstacle.                                          |
|                                                          | Function disabled on Attendant Control Panel                            | Enable disabled function.                                     |
|                                                          | Actuators have no power<br>Defective actuators<br>Defective accumulator | Check the mains connection.<br>Notify the service department. |
|                                                          | Mains Plug inserted incorrectly                                         | Insert the Mains Plug correctly.                              |
|                                                          | Faulty Power Source.                                                    | Notify the service department.                                |
|                                                          | Faulty Control Element.                                                 | Notify the service department.                                |
| Lowering Backrest from the upright position not possible | Obstacle under the Backrest or in the drive mechanism.                  | Remove the obstacle                                           |
|                                                          | CPR Release Handle is defective                                         | Notify the service department.                                |
| Adjusting Siderails not possible                         | Obstacle in the Siderail Release Mechanism                              | Remove the obstacle.                                          |
|                                                          | Siderail Release Mechanism is defective.                                | Notify the service department.                                |
| Faulty brakes                                            | Obstacle blocking brakes mechanically                                   | Remove the obstacle.                                          |
|                                                          | The brake mechanism is defective                                        | Notify the service department                                 |

## 19 Maintenance



### **WARNING!**

#### **Risk of injury when working on the bed!**

- ▶ Ensure that the bed is disconnected from the mains power prior to installation, putting into service, maintenance and deinstallation.
- ▶ Ensure that the castors are locked prior to installation, putting into service, maintenance and deinstallation.



### **WARNING!**

#### **Risk of injury due to defective bed!**

- ▶ Have a defective bed repaired immediately.
- ▶ If the defect cannot be repaired, do not use the bed.



### **CAUTION!**

#### **Material damage due to incorrect maintenance!**

- ▶ Ensure that maintenance is performed exclusively by manufacturer's customer service or by authorised service personnel certified by the manufacturer.
- ▶ If the defect cannot be repaired, do not use the bed.

LINET ® recommends attaching the maintenance plaque to the bed.

### 19.1 Regular maintenance

- ▶ Check regularly movable parts for wear.
- ▶ Perform regularly visual check of the product (with delivery note if necessary).
- ▶ Ask service department of the manufacturer for addition of the original spare parts if some product parts are missing.
- ▶ Ask service department of the manufacturer for replacement of any damaged product parts by the original spare parts.
- ▶ Check that the accumulator is working properly. Disconnect the bed from the mains power to check signalisation of accumulator indicator according to the instructions for use.
- ▶ Have the accumulator replaced if it is not working properly.
- ▶ Check regularly that all accessories are working properly.
- ▶ Replace damaged accessories immediately.

### 19.2 Spare Parts

The serial label is located on the frame of the mattress support platform. The serial label contains information for claims and ordering replacement parts.

#### **Information about spare parts is available from:**

- Manufacturer's customer service
- Sales department

### 19.3 Safety Technical Checks



### **WARNING!**

#### **Risk of injury due to incorrect safety technical checks!**

- ▶ Ensure that safety technical checks are performed exclusively by manufacturer's customer service or by authorised service personnel certified by the manufacturer.
- ▶ Ensure that the safety technical checks are recorded in the service and maintenance log.

**Safety technical check of the medical bed must be performed at least once every 12 months.**

The procedure for performing the safety technical check is stipulated in EN 62353:2014.

**NOTE** On request, the manufacturer will provide service documentation (e.g. circuit diagrams, component part lists, descriptions, calibration instructions etc.) for service personnel for the repair of ME equipment designated by the manufacturer as repairable by service personnel.

## 20 Disposal

### 20.1 Environment Protection

The company LINET® is aware of the importance of environmental protection for future generations. Within this company the environmental management system is applied in accordance with the internationally agreed standard ISO 14001. The compliance with this standard is annually tested by the external audit executed by an authorised company. Based on the Directive No. 2002/96/EC (Directive **WEEE** - Waste, Electric and Electronic Equipments) the company LINET, s. r. o. is registered in the List of Electric and Electronic Equipment Producers (**Seznam výrobců elektrozařízení**) on the Ministry of the Environment of the Czech Republic (Ministerstvo životního prostředí).

Materials used in this product are not environmentally hazardous. LINET® products meet valid requirements of national and European legislation in the areas of **RoHS** and **REACH**, so they do not contain any prohibited substances in excess quantities. None of the wooden parts is made of tropical wood (such as mahogany, rosewood, ebony, teak etc.) or made of timber from the Amazon region or from similar rainforests. Product noise (sound pressure level) meets requirements of the regulations for the protection of public health against undesirable effects of noise and vibration in protected interior spaces of buildings (according to standard IEC 60601-2-52). Used packaging materials are in accordance with requirements of the Packaging Act (**Zákon o obalech**). For disposal of packaging materials after installation of products contact your sales representative or manufacturer's customer service about the possibility of a free take-back of packaging through an authorized company (more details on [www.linnet.cz](http://www.linnet.cz)).

### 20.2 Disposal

The main objective of the obligations arising from the European Directive No. 2012/19/EU on Waste, Electric and Electronic Equipments (nationally regulated in Act No. 185/2001 Coll. as amended. On Waste and in Decree of the Ministry of the Environment No. 352/2005 Coll. as amended), is to increase the re-use, material recovery and recovery of electric and electronic equipment at the required level, thereby avoiding the production of waste and thereby avoiding the possible harmful effects of hazardous substances contained in electric and electronic equipment on human health and the environment. LINET® electric and electronic equipments that have a built-in battery or accumulator are designed so that the used batteries or accumulators can be safely removed by LINET® qualified service technicians. There is an information about its type on the built-in battery or accumulator.

#### 20.2.1 Within Europe

**To dispose of the electric and electronic equipment:**

- ▶ The electric and electronic equipment must not be disposed of as household waste.
- ▶ Dispose of this equipment at designated collection points or take-back points.

**To dispose of the other equipment:**

- ▶ The equipment must not be disposed of as household waste.
- ▶ Dispose of this equipment at designated collection points or take-back points.

LINET® participates in a collective system with take-back company REMA System (see [www.remasystem.cz/sberna-mista/](http://www.remasystem.cz/sberna-mista/)). By bringing electric and electronic equipment to a take-back point, you participate in recycling and you save primary raw material resources while protecting your environment from effects of unprofessional disposal.

#### 20.2.2 Outside Europe

- ▶ Dispose of the product or its components in accordance with local laws and regulations!
- ▶ Hire an approved waste disposal company for disposal!

---

## 21 Warranty

LINET® will only be held responsible for the safety and reliability of products that are regularly serviced, maintained and used in accordance with the safety guidelines.

**Should a serious defect arise that cannot be repaired during maintenance:**

- ▶ Do not continue to use the bed.

This product is covered by a 24-month warranty from the date of purchase. The warranty covers all material and manufacturing-related failures and errors. Failures and errors caused by incorrect use and external effects are not covered. Justified complaints will be fixed free of charge during the warranty period. Proof of purchase, with the date of purchase, is required for all warranty service. Our standard terms and conditions apply.

## 22 Standards and Regulations

Applied norms are stated on Declaration of Conformity.

The manufacturer adheres to a certified quality management system in compliance with the following standards:

- ISO 9001
- ISO 14001
- ISO 13485
- MDSAP (Medical Device Single Audit Program)